



LABWORKS Enterprise LIMS

User's Guide


Last update: Apr 28, 2021

Contents

Sign In	4
User Profile.....	5
Application Interface	6
Dashboard.....	7
Dashboard Selection.....	9
Default Dashboard	10
Dashboard Creation.....	11
Dashboard Copying	12
Dashboard Editing	13
Dashboard Tile Editing.....	15
Dashboard Sharing	17
Dashboard Deletion	19
Sample Management	20
Sample Login Groups.....	22
Create Samples	23
Sample Location Codes.....	26
Sample Special Info.....	27
Sample Analyses.....	28
Sample Attachments	30
Sample Report Recipients	31
Sample Coordinate Info.....	32
Work Queue Manager.....	33
Sidebar.....	36
Settings	37
Sample Explorer	39
Log Batch Explorer	41
Log Batch Details	42
QA/QC Batch Explorer.....	45
Collection Group Explorer	46
Explorer Settings.....	48
Explorer Filters.....	49
Quick Filters	50
Advanced Filters.....	51
Preset Filters	54
Views.....	55
Sample Details.....	56
Sample Editing	59
Sample Status History	61
Analysis Details.....	62
Analysis Editing	64
In-line Result Entry	66

Analysis Result History	67
Analysis Replicates	68
Enter Results	72
Edit / Review Samples	77
Receive Samples	78
Validate samples.....	80
Scan Samples.....	82
Audit Trail.....	83
Electronic signature	84
Analysis Legend	85
Analysis Types.....	85
Analysis States	85
Analysis Statuses.....	85
Analysis Result Colors	85
Collection Group Statuses	86
Checksum Violations	87

Sign In



Labworks Enterprise LIMS is a web application, so you need a web browser to access the system.

Supported desktop browsers:

- Google Chrome
- Safari

Supported mobile browsers:

- Chrome
- Safari

Please make sure that your browser is updated to the latest version.

The following information is required to access Labworks Enterprise LIMS:


1. URL-address of Labworks Enterprise LIMS in your network.
2. Database name, when there is more than one database available.
3. Your login and password in the selected database.

The first screen you see when you access the Labworks Enterprise LIMS for the first time is the **Sign In** screen. Once you fill up the form and click the "**Sign In**" button, the system will validate provided credentials and allow you to access the main interface of the application.

User Profile

In order to access various components of the Labworks Enterprise LIMS, a corresponding role must be assigned to your profile.

←



This user has not been assigned any licenses or permissions yet, please ask your administrator to set up this account. ×

User Profile Save

First name

Last name

Initials

Email

User groups

No value

User roles

No value

The system will display following notification in case if there are no roles assigned to your profile:


Please ask your administrator to assign certain roles to your profile.

Assigned roles are displayed on **User Profile** page under the "**User roles**" label.

Application Interface

Upon a successful sign in to the Labworks Enterprise LIMS, the user is navigated to the main area of the application interface. All screens in this area share the same layout consisting of two components:


1. **Application Header** is shared between all screens. It contains the global navigation menu and the user profile menu.
2. **Content Block** displays screen-specific information.




The global navigation menu is used to switch between three areas:

1. **Dashboard** area contains configurable tiles with a visual representation of various summary data.
2. **Sample Management** area provides various pages to manage samples: creating, viewing sample lists, searching, filtering, or inspecting particular samples.
3. **Configuration** area contains various system settings.

By clicking the **magnifier icon** you can open the **Search Form** and specify the text to search in the **Sample ID** field values of all available samples. Once the form is submitted, you are navigated to the **Sample Explorer** screen with a corresponding filter by **Sample ID** field applied, or the **Sample Details** screen, when there is only one sample found by this query.



The user profile menu is used to navigate to the **User Profile** screen and sign out of the application to end the session.




Dashboard

Dashboard is the first screen, the users with assigned roles see when they sign in to the system.


Every dashboard is a configurable set of tiles with a certain summary data.

The summary is generated only for the data you have access to, so **different users may see different data on the same dashboard**, depending on their access level.





The data displayed in the tiles updates automatically every five minutes.

You can force the update of a particular tile by clicking "**Refresh data**" option in the tile action menu.



To force the update of the whole dashboard at once, click the **refresh icon** next to the "**Last refresh**" label or the "**Refresh**" option in the dashboard action menu.




By clicking on charts inside the tiles you can navigate to the **Sample Explorer** screen with appropriate filters applied, to examine the data further.

Dashboard Selection

All dashboards can be divided into following groups by the level of access:

Icon	Access	Created by	Visible to	Configurable by
	Private (not shared)	Yourself	Yourself	Yourself
	Public (shared)	Yourself	All users	Yourself
	Public (shared)	Another user	All users	Dashboard creator

You can freely choose between your private dashboards or any of the public dashboards by clicking the dashboard title. The click propagates a window with a list of up to five recently viewed dashboards.



To access the full list of available dashboards, click "**Browse all dashboards**" button. The click opens the "**Dashboards**" sidebar, where you can search, filter and select any of available dashboards to be set as your dashboard.


Default Dashboard

The default is a dashboard displayed when you navigate to the **Dashboard** screen.

If you have no dashboard set as default, the global default dashboard will be displayed instead.


Since a public dashboard may be changed at any time by the owner, you can only set your own dashboard as a default one. To have the public dashboard set as your default, you have to create a private copy of that dashboard using the **Copy** action.

To set your own dashboard as a default, click the menu button on the right-hand side of the screen.




In the list of actions, find and click the "**Set as default**" option.

Now the dashboard is set as default, which is indicated by the check mark icon next to the dashboard title.




Dashboard Creation

If you wish to have a custom dashboard you can either copy and edit an existing one or create your own dashboard from scratch by clicking "**Add dashboard**" button on the right-hand side of the **Dashboard** screen.



In the pop-up window, enter the desired dashboard name and click "**Create**" button.


You will see an empty dashboard in the editing mode with the tile gallery opened on the left-hand side of the screen.



Add and configure tiles to your preference and click "**Save**" button to complete dashboard creation process.


Dashboard Copying

To copy an existing dashboard, select it in the list of dashboards, and click "**Copy**" option in the dashboard action menu.



In the pop-up window, enter the desired name for the copied dashboard and click "**Copy**" button.

You will see a copy of the dashboard in the editing mode with the tile gallery opened on the left-hand side of the screen.




Add and configure tiles to your preference and click "**Save**" button to complete dashboard copying process.

Dashboard Editing

You can edit only your own dashboards. If you wish to use an adjusted version of an existing dashboard of another user, you have to create a copy of that dashboard using the **Copy** action and edit the duplicate instead.

To edit a dashboard, select it in the list of dashboards, and click **"Edit"** option in the dashboard action menu. You will see the dashboard in the editing mode with the background filled with tile slots and the tile gallery opened on the left-hand side of the screen.



The tile gallery contains a grouped list of all available tile types with options to search through the tile names and collapse or expand a particular group.

While the dashboard is in the editing mode, you can perform following actions:

- **Rename the dashboard** in the text input field with the current dashboard name
- **Add tiles to the dashboard**
 - By selecting a tile in the tile gallery and clicking **"Add"** button that will add the tile to the next available space in the first row of the dashboard tile grid
 - By dragging a tile from the tile gallery with the left mouse button to the desired location in the dashboard

- **Relocate & reorder tiles** by dragging the tile within the dashboard bounds with the left mouse button
- **Resize tiles within the constraints of a particular tile type** by dragging tile border with the left mouse button

When you click "**Save**" button, the changes are applied to the dashboard and it leaves the editing mode.

The changes are reflected in the interface of other users once they refresh the dashboard, in case if the dashboard is public.

To discard all changes made to the dashboard, click "**Cancel**" button and confirm the action in the dialog window.


Dashboard Tile Editing

Tile editing and deletion can be performed from the tile action menu toggled by the button on the right-hand side of the tile.

To delete a tile, click the **"Delete"** option in the tile action menu of a particular tile.

To edit a tile, click the **"Configure"** option.

The users can edit only the tiles of their own dashboards. The **"Configure"** and **"Delete"** options are not available to everyone else. If you wish to use an adjusted version of an existing dashboard of another user, you have to create a copy of that dashboard using the **Copy** action, and edit the duplicate instead.



The **"Configure"** option brings up a tile configuration window. The window contains different settings depending on the tile type:

Location Code - the tile data is filtered according to the specified location

Sample Statuses - the tile data is filtered according to selected sample statuses

Departments - the tile data is filtered according to selected departments

Time periods - the tile data is filtered according to the specified time period

Miscellaneous - special configuration options, not related to the sample management, e.g., the time zone and hour format settings for the clock tile

Samples Completed by Location Code

Note: tiles will only display samples completed up to 24 hours ago.

Time periods

☐ Today

☒ 1 week

☐ 2 weeks

☐ 3 weeks

☐ Month

Location code

ABC

Cancel

Apply

Clock

Format

☒ 12 hour format

☐ 24 hour format

Time zone

(UTC+11:00) Solomon Is., New Caledonia

Cancel

Apply

The availability of tile configuration options depends on your access to particular locations and departments. When you copy a dashboard with unavailable data, a corresponding error message will be displayed in the tile configuration window

When you click "**Apply**" button, the changes are applied to the dashboard and it leaves the editing mode.

All changes made to the tiles are reflected in the interface of other users once they refresh the dashboard, in case if the dashboard is public.


To discard all changes made to the dashboard, click "**Cancel**" button and confirm the action in the dialog window.

Dashboard Sharing

Shared dashboard becomes publicly available to other users of the system.

To share a dashboard, select it in the list of dashboards, and click "**Share**" option in the dashboard action menu.

The "**Share**" option in the dashboard action menu is replaced by "**Unshare**" option once the dashboard is shared.




All other actions are still available to the dashboard owner.

Any change made to the public dashboard is reflected in the interface of other users viewing this dashboard.

The summary is generated only for the data you have access to, so **different users may see different data on the same dashboard**, depending on their access level.

The shared dashboard can be found by any user in the complete list of dashboards. The only options available to other users in the dashboard action menu are:

- **Refresh**
- **Copy**



The only option available to other users in the tile action menu is "**Refresh data**".

Samples by Status

8836 samples total



⋮

- Configure
- Refresh data
- Delete

With the "**Unshare**" action you can reset the dashboard back to private.


Private and deleted dashboards become unavailable to other users.

Dashboard Deletion

You can delete only your own dashboard, regardless of it being private or public.

Note that deleted dashboard becomes unavailable to other users in case if it was public (shared with other users).

To delete a dashboard, select it in the list of dashboards, click "**Delete**" option in the dashboard action menu and confirm the action in the dialog window.




Sample Management


Sample Management is an area of the system where users can access detailed information about the samples.


There are six options under the "Sample Management" section in the header:


1. Create samples
2. Work queue manager
3. Sample explorer
4. Log batch explorer
5. QA/QC batch explorer
6. Collection group explorer

There is a set of buttons you will see on some of the screens that perform certain actions with the samples or analyses:


 - Open **Receive samples** screen for the selected samples. The drop-down menu next to the button allows to select the template beforehand. By default, the last used template should be open. In the context menu this action is applied when you select the **Receive** menu item.

 - Open **Edit / review samples** screen for the selected samples. The drop-down menu next to the button allows to select the template beforehand. By default, the last used template should be open. In the context menu this action is applied when you select the **Edit / review** menu item.

 - Open **Enter results** screen for the selected samples or analyses. The drop-down menu next to the button allows to select the template beforehand. By default, the last used template should be open. In the context menu this action is applied when you select the **Enter results** menu item.

 - Prints standard labels for the selected samples. In the context menu this action is applied when you select the **Print -> Labels** menu item.

Some sample actions may be hidden under the action drop-down menu:



- Refresh – Refresh data for the selected samples.

- Set sample status – Change the sample status to the selected one. The set of all available statuses, the set of statuses available for manual transitions, the status order, and colors are defined by the workflow architect. The desktop **State Management** application is used for configuring workflow architect templates.
- Sample operations -> Reorder – Create samples based on the parent sample properties. The initial status and number of reordered samples are requested before reordering.
- Sample operations -> Delete – Delete the selected samples from the system.
- Print -> Custody tracking report – Print custody tracking report based on a selected template from the submenu.
- Change department – Change the sample's departments to the selected one from the **Change Department** window.
- Move to collection group – Move selected samples to a collection group that you select or create in the **Move to Collection Group** window.
- Move to log batch – Move selected samples to a log batch that you select in the **Move to Log Batch** window.

There are also two buttons that help to find and add samples to the system.




- Open **Scan Samples** window.



- Open **Create Samples** screen. The drop-down menu next to the button allows to select the template beforehand.

Some analyses actions may be hidden under the action drop-down menu:



- Refresh – Refresh data for the selected analyses.
- Assign -> Analyst – Assign an analyst to the selected analyses via the **Assign Analyst** window.
- Assign -> Instrument – Assign an instrument to the selected analyses via the **Assign Instrument** window.
- Assign -> Due date – Assign a due date to the selected analyses via the **Assign Due Date** window.
- Result history – Open the **Result history** window for the selected analysis.
- Scratch – Scratch the selected analyses.

Sample Login Groups

The list of login groups can be accessed by clicking the menu button on the right-hand side of the **Create Samples** page and selecting the “Manage login groups” option.

Manage Login Groups						
	Name	Description	Modified user	Modified date	Category	
✕	ap-test	No value	APHAM	08/12/2020 06:51 PM	No value	⋮
✕	kg_min	No value	KG	06/18/2020 12:28 PM	No value	⋮
✕	ap-05-20	No value	NP	05/20/2020 04:56 PM	No value	⋮
✕	56	No value	NP	05/15/2020 10:00 PM	No value	⋮
✕	mobile-2	No value	NP	05/15/2020 09:07 PM	No value	⋮
✕	mobile	No value	NP	05/15/2020 09:05 PM	No value	⋮
✕	DATE-TIME-CASE-1-23	No value	NP	05/14/2020 09:18 PM	No value	⋮
✕	DATE-TIME-CASE-1-22	No value	NP	05/14/2020 09:00 PM	No value	⋮
✕	DATE-TIME-CASE-1-2	No value	NP	05/14/2020 05:56 PM	No value	⋮
✕	DATE-TIME-CASE-1-1	No value	NP	05/14/2020 05:53 PM	No value	⋮
1 2 3 4 5 ... 20 21 → Items per page 10 1-10 of 201 results						
						Cancel Load group ▾

You can sort and filter the login groups to find the one you need, change the name, description or category of a particular group by selecting the “Edit” option in the group actions menu and delete login groups.

To use a particular login group for sample login, click on it to mark it as the selected one, and then click the “Load group” button.

To load only the location codes mentioned in the group, click on the corresponding option under the load group button drop-down menu.

Create Samples

Create Samples screen is intended to provide a convenient way to add new samples to the system.

Order	Location code	A_TEST	AGENCY	Analyses	Attachment	b_TEST	Save as login group
1	@2			Analyses	Attachments		
2	@3			Analyses	Attachments		100
3	@4			Analyses	Attachments		100

The sample data should be filled in a spreadsheet with rows corresponding to particular samples and columns – to their fields. The spreadsheet supports MS Office Excel style operations.

By default, the spreadsheet is blank and in order to add a new sample you need to specify a location code. You can do this either by typing it in the location code field, or by clicking the plus icon in the field to open the **Select Location Codes** window where you can browse all location codes with additional data and additional options to ease the process of sample addition.

The set of fields in the spreadsheet is determined by the selected **Template** and can be changed anytime. Note that if you try to change the template while the spreadsheet is not blank, all sample data will be discarded upon confirmation.

However, you can save the sample data for later by clicking the menu button on the right-hand side and selecting the **“Save as login group”** option.

Save as Login Group

Name

Start typing

The field is required.

Description

Category

Start typing

Cancel Save

The “**Manage login groups**” option will navigate you to the **Manage Login Groups** window where you can update, delete or load any of the available login group data to the spreadsheet, except the collection and submission dates (as they depend on the template settings). The data can be loaded even if the currently selected template is different from the template used to save login group.

The fields may contain not just numbers, dates and text information, but also references to other entities.

You can assign analyses that should be performed with a certain sample, add attachments or specify contacts and mail groups that should receive notifications about the samples.

Once you have added all required information about the samples and clicked the “Create samples” button, you will see the **Create Samples** window where you will be asked to enter the information related to the whole batch of samples (log batch).

Create Samples

Sample department	AQUASPLASH	Collection group	
Log batch	210414193	Group name	210414193
Log batch template ⓘ	LBTEMPLATENAME	Owner	Public
Log batch description		Description	
firstPrompt	Browse file		
secondPrompt	crystal-results-8d3e6b0c-2... Add attachment		

Cancel Create

The log batch may have a log batch template when the sample login template allows it. If it does, you can select a log batch template, and the fields related to that template will be listed under the template name.

The log batch may also be assigned to a collection group on the right-hand side of the window. You are free to specify the details of the collection group as well as the owner.


When you are done adding all required information about the log batch, you can click “Create”.

1 Sample Has Been Successfully Created

You can view the newly created sample(s) in Explorer using the "View in explorer" button, modify them using the "Edit / review samples" button or enter results by using the "Enter results" button.

View log batch

Close View in explorer Edit / review samples Enter results



All information is now added to the system and you can view the details of the samples and their log batch (if it was created with a template). You can also update the samples or proceed to the analysis result entry either with the template used last, or the template that you selected in the drop-down list next to the corresponding button.

Samples may support printable reports depending on the template used during sample creation. This is indicated by the Print drop-down menu in the top-right corner of the Create Samples Success window. Click on the corresponding report to download it as PDF file.

Sample Location Codes

Select Location Codes window is needed to find and select multiple location codes in order to populate the sample login spreadsheet.

Select Location Codes

View selected: 3

Clear selection

X

Location code	Count	Localione	A test	Agency	B test	Csr numb	Discharge no	Disposa
<input type="text"/>		<input type="text"/>						
<input type="checkbox"/> WGW	0	Walnut GWT Facility	No value	IRWD	No value	No value	No value	ROUTI
<input type="checkbox"/> WGWTEFFC	0	Walnut GWT Facili...	No value	IRWD	No value	No value	No value	ROUTI
<input type="checkbox"/> WGWTEFFG	0	Walnut GWT Facili...	No value	IRWD	No value	No value	No value	ROUTI
<input checked="" type="checkbox"/> WGW	3	Walnut GWT Facili...	No value	IRWD	No value	No value	No value	ROUTI
<input checked="" type="checkbox"/> WGW	2	Walnut GWT Facili...	No value	IRWD	No value	No value	No value	ROUTI
<input checked="" type="checkbox"/> WGW	1	Peters Canyon Ch...	No value	IRWD	No value	No value	No value	ROUTI
<input type="checkbox"/> WGW	0	Peters Canyon Ch...	No value	IRWD	No value	No value	No value	ROUTI
<input type="checkbox"/> WHATAU	0	Whatney - Auto Ce...	No value	IRWD	No value	No value	No value	ROUTI
<input type="checkbox"/> WHATAU_G	0	Whatney - Auto Ce...	No value	IRWD	No value	No value	No value	ROUTI
<input type="checkbox"/> WHATBU	0	Whatney - Burroug...	No value	IRWD	No value	No value	No value	ROUTI

1461-1470 of 1479 locations

1 2 ... 144 145 146 147 148

← →

Cancel

Add

The “**Count**” field corresponds to the number of samples you wish to add to the sample login spreadsheet with this location code.

You can use sorting and filters to find the particular codes and preview all codes that you’ve selected with the “**View selected**” button in the top-right corner of the window.

Once you have selected all location codes that you need and specified the number of samples you want to add with a particular code, click the “**Add**” button, and the samples will be added to the very end of the sample login spreadsheet.

Sample Special Info

Special Info window contains additional sample fields.



Special Information Form Fields

Text_1

Cancel Apply

All fields are optional and inherit from the location code by default.

Enter the corresponding sample data and click the **“Apply”** button to save it.

Sample Analyses

Analyses window lets you select the analyses to be performed with a particular sample. You can filter and sort the selected analyses or remove them from the list.

Analyses

Enter analysis code

+ Add analysis

	Analysis code	Analysis description	Department	Method refere...	Hold time	Price	Invoinc...	Repor...	Reporting unit	Special inform...	Work units
✖	!2	!2	KG2	No value	0	0			No value	No value	0
✖	\$ICP	Metals By ICP/MS	METALS	EPA 200.8/...	262080	13.17			ug/L	No value	0
✖	@4	@4	DSR	No value	0	0			No value	No value	0
✖	AG_ICP	Silver ICP	METALS	EPA 200.8	262080	0			ug/L	No value	0
✖	AS_ICP	Arsenic ICP	METALS	EPA 200.8	262080	0			ug/L	No value	0
✖	BE_ICP	Beryllium ICP	METALS	EPA 200.8	262080	0			ug/l	No value	0
✖	CD_ICP	Cadmium ICP	METALS	EPA 200.8	262080	0			ug/l	No value	0
✖	CR_ICP	Chromium ICP	METALS	EPA 200.8	262080	0			ug/L	No value	0
✖	CU_ICP	Copper ICP	METALS	EPA 200.8	262080	0			ug/L	No value	0
✖	FE_ICP	Iron	METALS	EPA 200.8 ...	262080	0			ug/L	No value	0

1 2

Items per page

10

1-10 of 19 results

✖ Delete all analyses

Cancel

Apply

By default, the list of analyses is empty or populated with the default analyses for the sample location and you have to type its name in the field to the left and select it from the suggested ones, or click the “**Add analysis**” button in the top-right corner to open the **Add Analyses** window.

On the **Add Analyses** window you can sort the analyses, filter them either by fields or by their Matrix filters, and select the analyses you wish to add to the sample.

Add Analyses

View selected: 0

Matrix filters

From

To

From


To

From

To

<

You can enter an additional analysis information straight away by clicking on the analysis code. You will see a reduced version of the **Analysis Details** screen where you can fill some of the analysis fields even before the sample is added to the system.


BE_ICP - Beryllium ICP

Cancel
Apply

No value
Status

No value
Validation state

No value
Validation user

No value
Validation date

No value
Sample ID

@0
WQ ID

No value
Group ID


No value
Container ID

GENERAL

CALCULATIONS

ATTACHMENTS

REPLICATES

Method ref... EPA 200.8 

Assigned i...

Due date

Price

Assigned a...

Comments

Special information form fields

Close


If you wish to specify replicates for the analyses with components right away, you can open the **“Replicates”** tab and set the number of expected replicates as well as the logic for analysis result calculation.

Once you added all required analyses, click the **“Apply”** button.

29

Sample Attachments

Samples may have files attached to them, and they are handled on the **Attachments** window.



The Attachments window displays a table with the following data:


	File name	Type	Description
	Bar (1).png		No value
	Bar (2).png		No value

Buttons: Cancel, Apply

The Sample may inherit some attachments from the location. In this case, you can't delete or update them.

To attach a new file to the sample, click the **Add attachment** button in the top-right corner of the window, and select the file by clicking the **Browse file** button. You may add a description of the file if you wish.

When done, click the **Add** button and the attachment will be added to the list.



The Add Attachment dialog box contains the following fields and buttons:

- File: No value
- Browse file button
- Description:
- Buttons: Cancel, Add

You can filter, sort delete attachments or update their description. Once you added all required attachments, click the **Apply** button.

Sample Report Recipients

If you want someone to receive reports about a particular sample, you can use the **Report Recipients** window.

Report Recipients											+ Add recipient + Add reference		✕
✕	↑ Reference	Doctype	Addrtype	Primary	Email	Fax number	Postal address	Addressee	Test string	Test			
✕	IDENTICAL	ident	ident	<input type="checkbox"/>	ident@mail.com	No value	No value	No value	No value	No value			
✕	IDENTICAL	ident	ident	<input type="checkbox"/>	ident@mail.com	No value	No value	No value	No value	No value			
✕	MAILLIST	mailmailmail	mailmail	<input type="checkbox"/>	useruseruserus...	No value	No value	No value	No value	No value			
✕	MAILLIST	mailmail	mail	<input type="checkbox"/>	admin@mail.c...	No value	No value	No value	No value	No value			
✕	MAILLIST	mail	mail	<input checked="" type="checkbox"/>	admin@mail.c...	No value	No value	No value	No value	No value			
✕ Delete all recipients											Cancel		Apply

The **Report Recipients** window is available when you click on the “**Report recipients**” cell during the sample login. It contains a list of all users that should receive the reports.

To add users, click the “**Add recipient**” button on the top-right and select those that should receive the report. To add a group of users, use the “**Add reference**” button.

The users can be filtered and sorted or deleted from the list.

Once you added all recipients, click the “**Apply**” button.

Sample Coordinate Info

Coordinate info is a window for geographical sample data input.

The screenshot shows a web application window titled "Coordinate Info". On the left, there is a sidebar with the heading "Sample collection location" and a green pin icon labeled "Set point". Below this are three input fields: "Address", "Latitude", and "Longitude". The "Address" field has a small location pin icon to its right. The main area of the window is a map of a region in Nebraska, showing towns like Arthur, Tryon, Stapleton, North Platte, and others, along with Lake McConaughy State Recreation. The map has a "Map" tab selected and a "Satellite" tab. At the bottom left of the map area is a "Launch Google Maps" link. At the bottom right are "Cancel" and "Apply" buttons. The map data is dated 8/2021.

The window is available when you click on the **Coordinate Info** field value during the sample creation.

You can only specify the **Sample collection location** by setting a location of a green pin on the map. Use the address field to find the needed location quicker.

Click the **"Apply"** button to save the coordinates.

Work Queue Manager

Work Queue Manager is intended to view queue samples or analyses to quick access to regular daily work.

The page requires an appropriate privilege to be assigned to your profile. Once you have this privilege assigned, the **Work Queue Manager** page becomes available.

The page has two tabs: Samples and Analyses.

The **Sample** tab displays a list of available samples that have not been completed or completed no more than a certain number of days ago.

Status	Sample ID	WQ ID	Submitted date	Collected date	Sample department	Violation	Analyses completed
Waiting for Report (68)							
Waiting for Report	AE48494	WGWTRCB	01/22/2021 08:09 PM	01/22/2021 12:00 AM	AQUASPLASH	No	2 of 2
Waiting for Report	AE48401	206	12/08/2020 07:13 PM	12/08/2020 09:00 AM	AQUASPLASH	No	1 of 1
Waiting for Report	AE47198	@I123	06/01/2020 08:05 PM	06/01/2020 12:00 AM	AQUASPLASH	No	No value
Waiting for Report	AE47195	AP-LOCCODE	06/01/2020 03:26 PM	06/01/2020 12:00 AM	AQUASPLASH	Yes	3 of 3
Waiting for Report	AE47194	AP-LOCCODE	06/01/2020 11:56 AM	06/01/2020 12:00 AM	AQUASPLASH	No	1 of 1
Waiting for Report	AE47187	@I123	06/01/2020 03:42 AM	06/01/2020 12:00 AM	AQUASPLASH	No	No value
Waiting for Report	AE47186	@I123	06/01/2020 03:42 AM	06/01/2020 12:00 AM	AQUASPLASH	No	No value
Waiting for Report	AE47180	AP-LOCCODE	05/29/2020 08:59 AM	05/29/2020 12:00 AM	AQUASPLASH	No	1 of 1

The **Analyses** tab displays a list of pending analyses.

Due status	Sample ID	Analysis code	Sample department	Status	Violation	QA/QC batch ID	Assigned analyst	Samp collect exp	Analysis complete date
Overdue (100)									
Overdue	AE48229	AP_LONG_CODE_ANALYSIS...	AQUASPLASH	Waiting for Analy...	No	No value	No value	No value	12/01/2020 08:22 PM
Overdue	AE48229	PH_S	AQUASPLASH	Waiting for Analy...	No	No value	No value	No value	12/01/2020 08:22 PM
Overdue	AE48194	!#\$%&-'*~()~[]	AQUASPLASH	Waiting for Analy...	No	No value	No value	No value	09/28/2020 03:09 PM
Overdue	AE48194	!#\$%&-'*~()~[]<>	AQUASPLASH	Waiting for Analy...	No	No value	No value	No value	09/28/2020 03:09 PM
Overdue	AE48194	!I	AQUASPLASH	Waiting for Analy...	No	No value	No value	No value	09/28/2020 03:09 PM
Overdue	AE48194	!INDE_QQQQQQQQQQQQ...	AQUASPLASH	Waiting for Analy...	No	No value	No value	No value	09/28/2020 03:09 PM

Colored columns can be collapsed/expanded.

The data in both tables can be sorted by clicking the field labels in the headers of the tables. Subsequent clicks on the same fields reverse the sorting direction.

Default order:

- **Samples**
 - sample status in the order defined by the workflow architect,
 - descending order by sample ID.
- **Analyses**
 - due status – first overdue, and then closest to overdue statuses,
 - descending order by sample ID,
 - ascending alphabetical order by analysis code.

The data in both tables can be filtered by top filters on the page.

Possible filters:


The table's columns are resizable. To change a column width, you can grab the separator in the table header and release it when the width suits you.

You can select one or more samples/analyses in the lists and perform actions on them. All general actions on samples and analyses are described in the **Sample Management** section above. To view the selected sample/analyses, turn on the “**Show selected**” toggle.


The total number of sample/analyses displayed on the screen is shown above the table.


The set of displayed fields can be changed in the **Settings** window for both tabs.

You can navigate to the details of a particular sample or analysis by clicking on corresponding fields within the table.


You can click the  button to refresh the data on the page or configure automatic refreshing in the **Settings** window. The last update time is displayed near the refresh button.

Sidebar


To view details for focused samples and analyses, you can enable the sidebar option. When clicking the  button, the sidebar will be open.



For samples, it displays analysis data with the ability to enter/edit analysis results, and status history data.




For analyses, it displays result history data.



Settings

The set of fields displayed in the **Work Queue Manager** tables can be configured in the **Settings** window.

To access the window, click the  button near the title.

The field configuration requires an appropriate privilege to be assigned to your profile. Once you have this privilege assigned, the **Settings** window becomes available in the edit mode.

Settings

SAMPLES

ANALYSES

GENERAL

+ Add attribute

Number	Field name
1.	Sample ID
<div>×</div> <div>2. ↓</div>	<div>Priority</div> <div>▼</div>
<div>×</div> <div>3. ↑ ↓</div>	<div>Collected date</div> <div>▼</div>
<div>×</div> <div>4. ↑ ↓</div>	<div>Sample department</div> <div>▼</div>
<div>×</div> <div>5. ↑ ↓</div>	<div>Analyses completed</div> <div>▼</div>
<div>×</div> <div>6. ↑</div>	<div>Violation</div> <div>▼</div>

Cancel

Save

On the **Samples** and **Analyses** tabs, you can change the order of the fields, remove certain fields from the appropriate table or add new ones.

On the **General** tab, you can configure auto-refresh settings. If auto-refresh is turned on, the **Work Queue Manager** page is automatically refreshed based on the specified refresh period.

Settings

SAMPLES

ANALYSES

GENERAL

Auto refresh

5 min

10 min

15 min

30 min

45 min

1 hour

7

min

Cancel

Save

"**Cancel**" button discards all changes in the settings.

"**Save**" button applies the settings and reloads the table.

Other users are not affected by your settings.

Sample Explorer

Sample Explorer is one of the main screens of the **Sample Management** area.

All available samples are displayed in rows with options to sort and filter the samples according to certain criteria.

The data in the table can be sorted by clicking the field labels in the headers of the main or nested tables.

Subsequent clicks on the same fields reverse the sorting direction.


Default order:

- **Samples** - descending chronological order by submission date
- **Analyses** - serial number of the analysis in the sample

The set of displayed fields can be changed in the "**Settings**" window.

Please refer to explorer filters for details about filtering options.

You can navigate to the details of a particular sample, analysis or analysis results by clicking on corresponding fields within the table.

LABWORKS
ENTERPRISE LIMS


DASHBOARD


SAMPLE MANAGEMENT

CONFIGURATION


Q


Labworks Fuser
Demo database






Sample Explorer  View

All



 Filters

 Sample ID	 Sample location	Requested date	Ok to collect	Status	Has replicates	Sample department	Priority
<input type="checkbox"/> AA82159	ZUMA 16080	No value	No value	Completed	No	No value	No value
<input type="checkbox"/> AA76231	Zone III P.S.	No value	No value	Completed	No	No value	No value
<input type="checkbox"/> AA77811	Zone D Pump Station Vault	No value	No value	Completed	No	No value	No value
<input type="checkbox"/> AC62929	Zone D - G Pump Station	03/25/2010 11:12 PM	Yes	Completed	No	NONPOTABLE	Medium
<input type="checkbox"/> AA24438	Zone B BPSfiber Influent	No value	No value	Completed	No	No value	No value
<input type="checkbox"/> AA24439	Zone B BPSfiber Effluent	No value	No value	Completed	No	No value	No value
<input type="checkbox"/> AD89986	Zone A/B Pump Station	06/03/2015 12:10 PM	Yes	Completed	No	NONPOTABLE	Low
<input type="checkbox"/> AD91191	Zone A/B Pump Station	06/17/2015 12:10 PM	Yes	Completed	No	NONPOTABLE	Low
<input type="checkbox"/> AD46993	Zone A/B Pump Station	10/16/2013 01:10 PM	Yes	Completed	No	NONPOTABLE	Low
<input type="checkbox"/> AE02502	Zone A/B Pump Station	11/18/2015 12:10 PM	Yes	Completed	No	NONPOTABLE	Low
<input type="checkbox"/> AD52305	Zone A/B Pump Station	01/15/2014 01:10 PM	Yes	Completed	No	NONPOTABLE	Low
<input type="checkbox"/> AE19497	Zone A/B Pump Station	08/03/2016 12:10 PM	Yes	Completed	No	NONPOTABLE	Low
<input type="checkbox"/> AD68839	Zone A/B Pump Station	08/20/2014 01:10 PM	Yes	Completed	No	NONPOTABLE	Low
<input type="checkbox"/> AE13092	Zone A/B Pump Station	05/04/2016 12:10 PM	Yes	Completed	No	NONPOTABLE	Low
<input type="checkbox"/> AE23555	Zone A/B Pump Station	10/05/2016 12:10 PM	Yes	Completed	No	NONPOTABLE	Low
<input type="checkbox"/> AD43786	Zone A/B Pump Station	09/04/2013 01:10 PM	Yes	Completed	No	NONPOTABLE	Low
<input type="checkbox"/> AD95448	Zone A/B Pump Station	08/19/2015 12:10 PM	Yes	Completed	No	NONPOTABLE	Low
<input type="checkbox"/> AE04740	Zone A/B Pump Station	01/06/2016 12:10 PM	Yes	Completed	No	NONPOTABLE	Low
<input type="checkbox"/> AD56770	Zone A/B Pump Station	03/19/2014 01:10 PM	Yes	Completed	No	NONPOTABLE	Low
<input type="checkbox"/> AD44655	Zone A/B Pump Station	09/18/2013 01:10 PM	Yes	Completed	No	NONPOTABLE	Low

1 2 3 4 5 ... 16963 16964 

Items per page

25

1-25 of 424085 samples

Samples selected: 0

When you click on the sample row, the element expands revealing the list of analyses performed with the sample.

<input type="checkbox"/>	AA82159	ZUMA 16080	No value	No value	Completed	No	No value
	Analysis code	Analysis name	Result	Start date	End date		
	C24_TCOL	Total Coliform (Colilert)	Absent	08/27/2001 08:43 PM	08/27/2001 08:43 PM		
	EC	Electrical Conductivity	408	08/24/2001 10:55 PM	08/24/2001 10:55 PM		
	TURB	Turbidity	0.709	08/24/2001 10:16 PM	08/24/2001 10:16 PM		
	ODOR	Odor - Threshold Odor Test	1	08/24/2001 10:37 PM	08/24/2001 10:37 PM		
	ODORDESC	Odor Description	Pipe Soap	08/24/2001 10:37 PM	08/24/2001 10:37 PM		
1 2	→	Items per page	5	1-5 of 10 analyses			

You can change the result of a particular analysis in the in-line mode according to the Inline result entry procedure or select a number of samples and fill their data on the enter result table by clicking the **"Enter Results"** button.

The **Analysis Replicates** screen of the single-component analysis can be accessed by clicking the value under the Replicate data.

The click on this value for multi-component analysis will open a window with a list of analysis components and their results. You can navigate to a particular component or browse its replicates from this window.

\$5252
✕

Component name	Result	
Benzy butyl phthalate	14	
Bis(2-ethylhexyl)phthalate	41	
Di-n-octyl phthalate	20	

Cancel
Save


You can select one or more samples in the list and perform actions on them. All general actions on samples are described in the **Sample Management** section above.

Log Batch Explorer

Log Batch Explorer has functionality similar to the sample explorer with exception that the samples are grouped in log batches and there are additional ways to sort and filter the table by the log batch fields.

LABWORKS		ENTERPRISE LIMS		DASHBOARD		SAMPLE MANAGEMENT		CONFIGURATION		Labworks Fuser		Demo database	
Log Batch Explorer										View All			
Filters													
Batch name		Last modified date		Total count of samples									
Unbatched		01/14/2021 07:20 PM		369									
P120709		04/13/2021 05:55 PM		236									
P101308		03/10/2011 08:03 PM		214									
P111008		03/10/2011 08:03 PM		208									
P090208		04/12/2021 10:56 AM		205									
P100713		12/22/2017 07:26 PM		203									
P080910		12/22/2017 06:23 PM		201									
P060908		03/10/2011 08:03 PM		199									
P040714		12/22/2017 07:38 PM		197									
P100812		12/16/2014 07:43 PM		196									
P100614		12/16/2017 05:12 PM		196									
P070813		12/22/2017 07:26 PM		194									
P110810		12/22/2017 06:23 PM		193									
P081009		06/15/2011 07:04 PM		193									
P080210		02/24/2016 10:55 PM		192									
P020711		12/22/2017 06:23 PM		192									
P022210		06/15/2011 07:04 PM		192									
P090808		03/10/2011 08:03 PM		191									
P102008		03/10/2011 08:03 PM		191									
P040813		12/22/2017 07:26 PM		191									
1 2 3 4 5 ... 1512 1513													
Items per page 25 1-25 of 37809 batches													
Samples selected: 0													

When you click on the batch row, the element expands revealing the list of samples in that batch. The expanded element works the same way as the list of samples on the **Sample Explorer** screen.

P120709		04/13/2021 05:55 PM		236			
	<input type="checkbox"/> Sample ID	Sample location	Requested date	Ok to collect	Status	Has replicates	Sample department
<input type="checkbox"/>	AC56267	Benner Reservoir-R43	12/13/2009 08:42 PM	Yes	Completed	No	POTABLE
<input type="checkbox"/>	AC56247	Benner Reservoir R43	12/12/2009 03:46 PM	Yes	Completed	No	POTABLE
<input type="checkbox"/>	AC56231	R43 Benner Reservoir	12/11/2009 06:37 PM	Yes	Completed	No	POTABLE
<input type="checkbox"/>	AC56185	Well #76, Yale @ I5 (IDP/PT...	12/10/2009 09:48 PM	Yes	Completed	No	POTABLE
<input type="checkbox"/>	AC56186	Well #77, Irvine High School...	12/10/2009 09:48 PM	Yes	Completed	No	POTABLE
<input type="checkbox"/>	AC56188	Well #107, Culver @ Trabuc...	12/10/2009 09:48 PM	Yes	Completed	No	POTABLE
<input type="checkbox"/>	AC56187	Well #115, Commerce (IDP/...	12/10/2009 09:48 PM	Yes	Completed	No	POTABLE
<input type="checkbox"/>	AC56184	Well 110, Walnut @ I5 (IDP/...	12/10/2009 09:48 PM	Yes	Completed	No	POTABLE
<input type="checkbox"/>	AC56189	Well Head at OPA Well 3	12/10/2009 09:48 PM	Yes	Completed	No	POTABLE

You can select one or more samples in the list and perform actions on them. All general actions on samples are described in the **Sample Management** section above.

The actions can be applied both to individually selected samples and to the entire group via the “dot-dot-dot” menu.

Some log batches may contain various information depending on the template it was created with. This is indicated by an underline when you hover over the Batch name. Click on the batch name to open the **Log Batch Details** window.

Log Batch Details

Certain log batch templates allow creation of complex log batches that may contain additional information related to all samples in the batch. The **Log Batch Details** window contains all information related to a particular log batch under the tabs.


The **“General”** tab contains the main log batch information along with the special forms data.

There are two types of special information forms:

1. Specified according to the log batch template with preset labels (prompts)
2. Unspecified special forms

Users can leave comments for the batch. The comments will appear in the bottom of the **“General”** tab.

200303070502 - Trafalgar square fountain water samples Edit

 APHAM
Modified user

03/03/2020 03:12 PM
Modified date

GENERAL

SAMPLES

PARAMETERS

ATTACHMENTS

REPORT RECIPIENTS

Log batch 200303070502

Log batch description Trafalgar square fountain water samples

▼ Special information forms

+ Add special information form

▶ Swamp chem

▼ TEST

Observed Color No value

Observed Odor No value

Observed Turbidity ... No value

Dewatering Activity ... 1

Receiving Water 3

Observed Suspen... . No value

Observed Oil and 2


▼ Comments

No value

Close

The **“Samples”** tab displays a list of samples of that log batch with sorting and filtering options.

200303070502 - Trafalgar square fountain water samples Edit

 APHAM
Modified user

03/03/2020 03:12 PM
Modified date

GENERAL

SAMPLES

PARAMETERS

ATTACHMENTS

REPORT RECIPIENTS

Sample ID

Group ID

WQ ID

Collected date

Submitted date

Modified date

Login user

Status

AE36409

2003030...

AT1_AN1

03/03/2020 03...

03/03/2020 03...

03/03/2020 03...

APHAM

Waitin...


Items per page 10

1-1 of 1 samples

Close

The **“Parameters”** tab contains custom fields of for this particular log batch. You can add, remove parameters, and set labels for the unspecified parameters.

200303070502 - Trafalgar square fountain water samples Edit

 APHAM
Modified user

03/03/2020 03:12 PM
Modified date

GENERAL

SAMPLES

PARAMETERS

ATTACHMENTS

REPORT RECIPIENTS

+ Add parameter

TEST

Active

Expected d...

03/05/2020 03:09 PM

Close

The “**Attachments**” tab contains all attachments related to the log batch in two groups:

3. Attachments specified according to the log batch template
4. Unspecified attachments

GENERAL	SAMPLES	PARAMETERS	ATTACHMENTS	REPORT RECIPIENTS
Specified + Add attachment				
<input type="checkbox"/>	<input type="text" value="Prompt"/>	<input type="text" value="File name"/>	<input type="text" value="Type"/>	<input type="text" value="Modified date"/>
	the first doc	1582243197423.jpg	Copy	03/03/2020 03:11 PM
Unspecified + Add attachment				
<input type="checkbox"/>	<input type="text" value="Prompt"/>	<input type="text" value="File name"/>	<input type="text" value="Type"/>	<input type="text" value="Modified date"/>
	No value	book.pdf	Copy	03/03/2020 03:11 PM

The “**Report recipients**” tab is similar to the **Sample Report Recipients** window, as it contains a list of all users that should receive the reports about the log batch.

GENERAL	SAMPLES	PARAMETERS	ATTACHMENTS	REPORT RECIPIENTS
+ Add recipient + Add reference				
<input type="checkbox"/>	<input type="text" value="Reference"/>	<input type="text" value="Doctype"/>	<input type="text" value="Addrtype"/>	<input type="text" value="Primary"/>
	IDENTICAL	ident	ident	<input type="checkbox"/>

You can update the information displayed in the tabs or leave a comment by entering the “Edit” mode. To do this, you can click the “**Edit**” button at the top-right corner of the screen.

The “Edit” mode enables automatically when you update the batch by clicking on the following buttons:

- **Add special information form** button on the **General** tab,
- **Add parameter** button on the **Parameters** tab,
- **Add attachment** button on the **Attachments** tab,
- **Add recipient / Add reference** buttons on the **Report recipients** tab;

The “Edit” mode is indicated by the “**Save**” and “**Cancel**” button at the top-right corner of the window as well as some fields becoming editable.

200303070502 - Trafalgar square fountain water samples

Cancel
Save

APHAM

Modified user

03/03/2020 03:12 PM

Modified date

GENERAL

SAMPLES

PARAMETERS

ATTACHMENTS

REPORT RECIPIENTS

Log batch

200303070502

Log batch description

Trafalgar square fountain water samples

Special information forms

+ Add special information form

Swamp chem

TEST

Observed C...

Observed O...

Observed S...

Observed T...

Dewatering...

Observed O...

Receiving ...

3

1

2

Comments

APHAM

First comment

Your comments

Add

Close


Click the **“Save”** button to apply the changes, and **“Cancel”** – to revert them.

All changes made to the log batch are reflected in the **Audit Trail**, available under the drop-down menu on the top-right.

QA/QC Batch Explorer

QA/QC Batch Explorer inherits the functions of log batch explorer, but displays the samples in QA/QC batches instead of log batches.

Note that the same sample may be a member of multiple QA/QC batches and represented in the table multiple times. These samples work synchronously when you perform certain actions with them, such as inline result entry.

LABWORKS
ENTERPRISE LIMS


DASHBOARD


SAMPLE MANAGEMENT

CONFIGURATION



Q


Labworks Fuser
Demo database





QA/QC Batch Explorer  View

All

 Filters

 Batch name	↓ Last modified date	Total count of samples	Batch number	Batch test	Batch date
Unbatched	04/15/2021 08:59 PM	67888	No value	No value	No value
TSS-181828	04/15/2021 09:59 AM	7	181828	TSS	05/02/2016 04:0
TSS-145645	04/07/2021 01:46 PM	10	145645	TSS	04/19/2013 08:5
PH-145636	04/07/2021 01:46 PM	3	145636	PH	04/19/2013 07:0
VSS-145647	04/07/2021 01:46 PM	10	145647	VSS	04/19/2013 08:5
ALKDA-145648	04/07/2021 01:46 PM	5	145648	ALKDA	04/19/2013 09:0
TURB-185463	04/02/2021 11:43 AM	10	185463	TURB	09/07/2016 10:1
NO2DA-185466	04/02/2021 11:43 AM	7	185466	NO2DA	09/07/2016 10:2
EC-185484	04/02/2021 11:43 AM	10	185484	EC	09/08/2016 06:1
AP-ANALYSIS-190242	01/14/2021 07:20 PM	10	190242	AP-ANALYSIS	02/28/2020 10:5
14DIOX-185475	01/14/2021 06:16 PM	4	185475	14DIOX	09/08/2016 03:5
ALK-2879	12/23/2020 04:32 PM	5	2879	ALK	11/03/1997 11:3
NH3_E-2865	12/23/2020 04:32 PM	8	2865	NH3_E	11/03/1997 07:3
BOD-3030	12/23/2020 04:32 PM	3	3030	BOD	11/11/1997 05:5
CODD-2870	12/23/2020 04:32 PM	2	2870	CODD	11/03/1997 11:1
TSS-2866	12/23/2020 04:32 PM	9	2866	TSS	11/03/1997 07:3
VSS-2867	12/23/2020 04:32 PM	7	2867	VSS	11/03/1997 07:3
CODT-2871	12/23/2020 04:32 PM	8	2871	CODT	11/03/1997 11:2
IBACTI-190177	12/10/2020 09:07 AM	30	190177	IBACTI	02/01/2018 09:5
AP_LONG_CODE_ANALYSIS_QQ-190165	12/10/2020 09:07 AM	5	190165	AP_LONG_CODE_ANALYSIS_QQ	01/26/2018 10:3

1 2 3 4 5 ... 7006 7007  Items per page

25

 1-25 of 175170 batches Samples selected: 0

You can select one or more samples in the list and perform actions on them. All general actions on samples are described in the **Sample Management** section above.

The actions can be applied both to individually selected samples and to the entire group via the “dot-dot-dot” menu.

Collection Group Explorer

Collection group explorer displays the information about the samples in a similar way to the **Sample Explorer**, except that all samples are shown under their respective collection groups. There are also special options to filter and sort by collection group fields.

To view the samples of a particular collection group, click on it to expand the list of samples. The list of samples supports sorting and filtering just like on the **Sample Explorer** page.

Collection Group Explorer

Filters: Sample created date: Last 30 days

Name	Description	Owner	Created date	Modified date	Ready to Collect
210415200	No value	Public	04/15/2021 08:59 PM	No value	1

Sample ID: AE48513, Sample location: No value, Requested date: 04/15/2021 08:58 PM, Ok to collect: Yes, Status: Ready to Coll..., Has replicates: No, Sample department: AQUASPLASH, Priority: High

Items per page: 25, 1-1 of 1 samples

Sample ID	Sample location	Requested date	Ok to collect	Status	Has replicates	Sample department	Priority
210415199	No value	04/15/2021 08:58 PM	Yes	Ready to Coll...	No	AQUASPLASH	High
210415195	No value	04/15/2021 08:51 PM	Yes	Ready to Coll...	No	AQUASPLASH	High
210414194	No value	04/14/2021 03:19 PM	Yes	Ready to Coll...	No	AQUASPLASH	High
210407190	No value	04/07/2021 05:38 PM	Yes	Ready to Coll...	No	AQUASPLASH	High
210407188	66666666666666666666	04/07/2021 05:35 PM	Yes	Ready to Coll...	No	AQUASPLASH	High
210401186	No value	04/01/2021 05:04 PM	Yes	Ready to Coll...	No	AQUASPLASH	High

Items per page: 25, 1-7 of 7 groups, Samples selected: 0

There is a “dot-dot-dot” menu on the right-hand side of every collection group that contains the list of actions available to the user.

Name	Description	Owner	Created date	Modified date	Ready to Collect	Collected	Not Collected	Status
120220009	No value	Public	12/02/2020 03:19 PM	No value	1	0	0	Audit trail
120120007	No value	Public	12/01/2020 06:21 PM	No value	1	0	0	Edit group
112520155	No value	Public	11/25/2020 05:53 PM	No value	0	0	0	Ungroup
112520153	No value	Public	11/25/2020 05:03 PM	No value	2	0	0	Edit/Review
112520152	No value	Public	11/25/2020 05:02 PM	No value	1	0	0	Receive
112520143	No value	Public	11/25/2020 04:26 PM	11/25/2020 05:04 PM	0	0	0	Set sample state
112520142	No value	Public	11/25/2020 04:25 PM	No value	0	0	0	Sample operations
112520138	No value	AA-TEST	11/25/2020 11:47 AM	No value	1	0	0	Print
112420137	No value	JR	11/24/2020 10:32 AM	No value	0	0	0	Enter results
111820128	No value	Public	11/18/2020 05:22 PM	No value	0	0	1	Change department
111320122	No value	JR	11/13/2020 01:51 PM	No value	0	0	1	Move to collection group
111320121_jr	No value	JR	11/13/2020 01:43 PM	No value	0	0	0	Move to log batch
111320119_test_audit_trail	No value	JR	11/13/2020 01:07 PM	No value	0	0	1	Mark received

Apart from the regular actions, there are also ones available to collection groups only:

- **Edit group** opens a window where you can change the collection group attributes such as name, description and owner.
- **Audit trail** opens the **Audit Trail** for the particular collection group.
- **Ungroup** action deletes the collection group removing all samples from it as if they were added without a collection group. However, the sample data is not lost, as all samples of the same log batch without a group are shown in the list under a fake collection group that has no data and called “Ungrouped”.

You can select one or more samples in the list and perform actions on them. All general actions on samples are described in the **Sample Management** section above.

The actions can be applied both to individually selected samples and to the entire group via the “dot-dot-dot” menu.

Explorer Settings

The set of fields displayed in the **Sample Management Explorer** table can be configured in the "**Settings**" window.

To access the window, click the  **Column options** button above the sample explorer table.

The field configuration requires an appropriate privilege to be assigned to your profile. Once you have this privilege assigned, the "**Settings**" window becomes available in the edit mode.

Settings

+ Add field

X

Order	Field ID	Field name
1.	Sample ID	Sample ID
<div>X</div> 2. <div>↓</div>	Sample location	<div>Sample location</div>
<div>X</div> 3. <div>↑ ↓</div>	Requested date	<div>Requested date</div>
<div>X</div> 4. <div>↑ ↓</div>	OK to collect	<div>Ok to collect</div>
<div>X</div> 5. <div>↑ ↓</div>	Status	<div>Status</div>
<div>X</div> 6. <div>↑ ↓</div>	Has replicates	<div>Has replicates</div>
<div>X</div> 7. <div>↑ ↓</div>	Sample department	<div>Sample department</div>
<div>X</div> 8. <div>↑ ↓</div>	Priority	<div>Priority</div>
<div>X</div> 9. <div>↑ ↓</div>	Priority	<div>Priority</div>
<div>X</div> 10. <div>↑ ↓</div>	Group ID	<div>Group ID</div>
<div>X</div> 11. <div>↑</div>	Price list	<div>Price list</div>

Cancel

Save

Here you can change the order of the fields, rename them, remove certain fields from the table or add new ones.

"**Cancel**" button discards all changes in the settings.

"**Save**" button applies the settings and reloads the table.

Other users are not affected by your settings.

Explorer Filters

The data in the **Explorer** area is displayed according to the filters.

You can narrow down the displayed data using four independent types of filters:

- **Quick Filters** available in the quick filter bar opened by clicking the **gray funnel icon** in the table header.
- **Views** available under the **View** drop-down list.
- **Advanced Filters** available in the sidebar opened by clicking the **white funnel icon** above the table.
- **Preset Filters** available in the very bottom of the sidebar.

Quick Filters

Quick Filters allow you to apply simple filters directly from the tables.

To filter the elements, expand the quick filter bar by clicking the **funnel icon** on the left-hand side of the table header, specify the filters, and apply them by clicking **check mark icon** in the quick filter bar. Moreover, the filters will be applied as soon as you move the focus to another control or just press Enter.

The **cross icon** in the quick filter bar clears all quick filters.


 <input type="checkbox"/> Sample ID	Sample location	Requested date	Ok to collect	Status	 Has replicates	Sample department
  <input type="text"/>	<input type="text"/>	<input type="text"/> 	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

Advanced Filters


Advanced filters give you extended control over the table by allowing to filter the data by hidden fields.

The **Advanced Filters** are configured in the sidebar, where all available filters are listed in expandable groups.

The sidebar is toggled by the **funnel icon** above the table.



By clicking the **pencil icon** next to the group name, you can select filters that should appear in the sidebar.

A dialog box titled "Sample Field Filters Configuration" with a close button (X) in the top right corner. It contains a list of filter fields with checkboxes. The first field, "Bac-T-number", is selected with an orange square. The other fields are "Bill rate", "Billing Address line 1", "Billing Address line 2", "Billing Address line 3", "Billing Address line 4", "Billing Address line 5", "Comments mask", and "CTReport", all with empty checkboxes. At the bottom, there is a "Reset" button with a red X icon, a "Cancel" button, and a "Save" button.

Sample Field Filters Configuration

- ☒ Bac-T-number
- ☐ Bill rate
- ☐ Billing Address line 1
- ☐ Billing Address line 2
- ☐ Billing Address line 3
- ☐ Billing Address line 4
- ☐ Billing Address line 5
- ☐ Comments mask
- ☐ CTReport

× Reset Cancel Save

The advanced filters are divided into four groups:

- **Sample filters** applied to the lists of samples
- **Analysis filters** applied to the lists of analyses
- **QA/QC batch filters** applied to the list of QA/QC batches (Appears on QA/QC batch explorer only)
- **Collection group filters** applied to the list of Collection groups (Appears on Collection group explorer only)

Sample and analysis filters are also divided into groups by field data type:

- **Date filters** applied to all fields of date-time data type
- **Field filters** applied to the fields of other data types

Expiration Date

Calendar day

☒ Today
☐ Number of last days

☐ Yesterday

☐ Last 7 days

☐ Last 15 days

☐ Last 30 days

☐ Custom period

12/09/2020

12/09/2020

Calendar week

☐ This week
 12/09/2020

☐ Last week

Calendar quarter

☐ This quarter

☐ Last quarter

Calendar year

☐ This year

☐ Last year

Cancel

Apply

Status

☐ Sample Recorded
☐ Not Ready

☐ Ready to Collect
☐ Collected

☐ Not Collected
☐ In Transit

☐ Received
☐ Waiting for Analyses

☐ Analyses completed
☐ Waiting for Validation

☐ Validation Completed
☐ Waiting for Report

☐ Report Completed
☐ Report Rejected

☐ Waiting for Invoice
☐ Invoice Completed

☐ Invoice Rejected
☐ Completed

☐ Retired
☐ Deleted

Cancel

Apply

Once you specify an advanced filter, it appears in the new chip above the explorer table.

Filters

Submitted date

Last 7 days

Sample status

Not Ready X

Collected X

Completed X

X Clear filters

You can adjust the filter by clicking on a particular chip or remove it with the **cross icon** inside the chip.

"Clear filters" button on the right-hand side of the screen removes all custom filters applying only the default ones instead.

Preset Filters

The **Preset Filters** contain most popular sets of filters. You can find them in the sidebar toggled by the **funnel icon** above the table.

Preset Filters

☐ Analyses count by departments

☐ Samples by state

☒ Samples by state and sample department

☐ Samples completed by location code

☐ Samples out of specification in date range

☐ Samples waiting for approval by departments

☐ Today's work – queued vs completed

☐ Top 10 location codes

☐ Top 15 viewed

☐ Last 15 viewed

☐ Samples by state and analysis department

☐ Samples completed

☐ Samples currently out of specification

☐ Samples waiting for approval

☐ Today's work overview

☐ Top 10 analyses

☐ Top 10 location codes by completed samples

Cancel

Apply

Once you apply one of the presets, this will be indicated above the explorer table.

 Filters

Preset

Top 10 location codes X

Location codes

@0 +9

Date


Last 30 days

You can adjust the filter by clicking on the corresponding chips or by opening the sidebar and changing the preset filter values there.

The dashboard tiles apply one of these presets automatically when you navigate to the sample explorer by clicking on the dashboard tile chart.

Views

The **Views** work like your personal preset filters that can be customized and saved for later. Views are saved separately for each sample management table and available only to the user who created them.



The default view "**All**" displays data about all samples.

You can create a new view by applying desired filters (quick or advanced) to the table and selecting the "**Save as new**" option under the drop-down menu next to the **View** field.




Once you provide the name and click the "**Create**" button, the new view becomes available as one of the **View** field options and applied as the active one.

You can update the custom view while it is active by adjusting filters and clicking the "**Save**" button appeared next to the **View** field.

To rename custom view, select the "**Rename**" option under the drop-down menu next to the **View** field.


To delete it, select the "**Delete**" option.



Sample Details

Sample Details screen displays all information about the sample and all related analyses.

The “dot-dot-dot” menu provides users with the ability to perform actions on samples. All general actions on samples are described in the **Sample Management** section above.



The "**Audit Trail**" option allows you to browse all changes made to the sample on the **Audit Trail** window.

The "**Status history**" option opens the **Status History** window with the information on sample status changes.

The sample data is displayed on multiple tabs:

General tab contains the main sample data that is organized in expandable.

A screenshot of the 'General' tab in the Sample Details interface. The interface shows a sidebar on the left with sample information: 'Completed Status', 'SJRES_CB WQ ID', 'San Joaquin Reservoir, Site C Bottom Sample location', 'Bonnie Steward Sample Collector', '12/30/2008 10:39 PM Collected date', and '03/10/2010 11:37 PM Submitted date'. The main area displays a form with fields for 'BSTEWARD Login user', 'Modified date' (04/14/2021 02:53 PM), 'Sample due date' (01/20/2009 08:00 AM), 'Priority' (Low), 'Sample department' (DATS), 'Number of labels' (1), 'Customer' (No value), and 'Collection group' (No value). There are also checkboxes for 'Report queue', 'Do validation option value', 'Do invoicing option value', and 'Is received'. At the bottom, there are expandable sections for 'User defined fields', 'Special information form fields', and 'QA/QC batch membership'. A 'Close' button is in the bottom right corner.

The click on the sample **Group ID** name opens the **Log Batch Details** page. If the log batch details do not exist, the **Log Batch Details** page will be open with default values. The new details will be added to the system after saving.

The related analyses are displayed under the **Analyses** tab in a list with filtering options. You can navigate to a particular analysis by clicking the analysis code, or enter the result directly in the list.

AC34273

Status history

Edit

Completed

Status

SJRES_CB

WQ ID

San Joaquin Reservoir, Site C Bottom

Sample location

Bonnie Steward

Sample Collector

12/30/2008 10:39 PM

Collected date

03/10/2010 11:37 PM

Submitted date

GENERAL

ANALYSES

ATTACHMENTS

COORDINATE

REPORT RECIPIENTS

Analysis code

Analysis name

Result

Start date

End date

Replicate ...

Analyst

Status

DO_FLD

Dissolved Oxygen (fi...

10.7

12/30/2008 10:39 PM

12/30/2008 10:39 PM

No

RPEREZ

Manually ent...

EC_FLD

Electrical Conductivi...

1200

12/30/2008 10:39 PM

12/30/2008 10:39 PM

No

RPEREZ

Manually ent...

TEMP_FLD

Temperature (field)

17.93

12/30/2008 10:39 PM

12/30/2008 10:39 PM

No

RPEREZ

Manually ent...

Items per page

10

1-3 of 3 analyses

Close

AC34273
Status history Edit

GENERAL	ANALYSES	ATTACHMENTS	COORDINATE	REPORT RECIPIENTS
<p>Completed Status</p> <hr/> <p>SJRES_CB WQ ID</p> <hr/> <p>San Joaquin Reservoir, Site C Bottom Sample location</p> <hr/> <p> Bonnie Steward Sample Collector</p> <hr/> <p>12/30/2008 10:39 PM Collected date</p> <hr/> <p>03/10/2010 11:37 PM Submitted date</p>	<div style="position: absolute; top: 10px; right: 10px; background: white; padding: 5px; border: 1px solid #ccc;"> Sample collection location Address Albuquerque International Support-Abq (ABQ), Sunport Blvd, Albuquerque, NM, USA Latitude 35.0433333 Longitude -106.6129085 </div>			
	<div style="position: absolute; top: 10px; right: 10px; background: white; padding: 5px; border: 1px solid #ccc;"> Field collection location Address Albuquerque International Support-Abq (ABQ), Sunport Blvd, Albuquerque, NM, USA Latitude 35.0433999 Longitude -106.6129085 0.005 mi Distance </div>			

AC34273

Status history 📄 ⚙️ 🔍 Edit ✕


GENERAL	ANALYSES	ATTACHMENTS	COORDINATE	REPORT RECIPIENTS					
SJRES_CB WQ ID San Joaquin Reservoir, Site C Bottom Sample location Bonnie Steward Sample Collector 12/30/2008 10:39 PM Collected date 03/10/2010 11:37 PM Submitted date									
+ Add recipient + Add reference									
🔗 Reference	Ddoctype	Addrtype	Primary	Email	Fax number	Postal address	Addressee	Test string	Test
BOOLTEST	boolT	boolT	■	boolT@mail.c...	No value	No value	No value	No value	No value
IDENTICAL	ident	ident	■	ident@mail.c...	No value	No value	No value	No value	No value


Close

Sample Editing

Sample Editing mode is enabled when you press the **"Edit"** button on the **Sample Details** screen, by double-clicking any of the editable fields or with the **"Add analysis/attachment/recipient/reference"** button in the tabs.

In the edit mode, certain fields of the sample on all tabs become editable and the **"Edit"** button is replaced with **"Save"** and **"Cancel"** buttons that apply or discard the changes respectively.

Every change can be discarded separately by the **"Revert"** button  next to the updated field.




The screenshot shows the 'Sample Details' screen for sample AE41352 in 'Edit' mode. The left sidebar contains a list of fields: 'Collected State', 'APLATT Login user', 'No value Sample Collector', '\$TEST_LOCCODE_0077 Location code', 'No value Location description', and '111220076 Collection group'. The main area has tabs for 'GENERAL', 'ANALYSES', 'ATTACHMENTS', 'COORDINATE', and 'REPORT RECIPIENTS'. The 'GENERAL' tab is active, showing fields like 'Log batch' (111220076), 'Submitted date' (11/12/2020 10:41 PM), 'Collected date' (11/12/2020 10:00 AM), 'Sample due date' (11/12/2020 10:00 AM), 'Modified date' (12/03/2020 04:55 PM), and 'Comments'. There are also checkboxes for 'Report queue', 'Do validation option value', 'Do invoicing option value', and 'Is received'. At the bottom right, there is a 'Manage membership' link and a 'Close' button.

The Login user and Sample collector can be selected in the user selection window opened by the corresponding buttons next to these fields.

The available users are displayed in a paginated list with sorting and filtering options.

To set a user as the login user or the sample collector, select him in the list and click the **"Apply"** button.







The screenshot shows the 'Login User' selection window. It has a search bar with a filter icon and a dropdown for 'Initials'. Below the search bar is a table with two columns: 'First and last name' and 'Test'. The table lists two users: 'AA-TEST' and 'TPK'. At the bottom, there is a pagination control showing 'Items per page' (10) and '1-2 of 2 users'. There are 'Cancel' and 'Apply' buttons at the bottom right.

The batches listed under the **QA/QC batch membership** group can be selected in the corresponding window opened by the "Manage membership" button.

The batches are displayed in a paginated list with sorting and filtering options.

You can select a one or many batches, preview the selection by clicking the **"Batches selected: X"** option, clear the selection with **"Clear selection"** button or apply the selection by clicking the **"Apply"** button.

QA/QC Batch		Batches selected: 2	✖ Clear selection	✕
  Batch name	Total count of samples			
<input type="checkbox"/> ZN_ICP-99746	2			
<input type="checkbox"/> ZN_ICP-190235	2			
<input type="checkbox"/> ZN_ICP-190219	2			
<input checked="" type="checkbox"/> X_ZN_S-190192	3			
<input checked="" type="checkbox"/> X_SE_S-190191	4			
<input type="checkbox"/> X_PB_S-190190	4			
<input type="checkbox"/> X_NI_S-190189	3			
<input type="checkbox"/> X_MO_S-190188	3			
<input type="checkbox"/> X_MBAS-68701	3			
<input type="checkbox"/> X_HG_S-190187	2			
1 2 3 4 5 ... 17754 17755  Items per page <input type="text" value="10"/> 1-10 of 177546 batches				
		 		

QA/QC Batch		View all	✖ Clear selection	✕
  Batch name	Total count of samples			
<input checked="" type="checkbox"/> X_ZN_S-190192	3			
<input checked="" type="checkbox"/> X_SE_S-190191	4			
Items per page <input type="text" value="10"/> 1-2 of 2 batches				
		 		

Sample Status History


The **Status History** screen is where you can see how the status of the sample has been changing over time.

The history of statuses is displayed in the chart in chronological order and color-coded according to the legend.

The set of sample statuses and their colors are defined by the workflow architect template. It is configured by the desktop State Management application.

The "**Completed**" status, defined by the default workflow architect template, is excluded from the chart when the sample has it as the final one.

You can see the exact time when the sample entered and left a certain status as well as its time period by hovering the cursor over a certain status in the chart.




You can see the audit trail where logs are filtered only by status changes.

Analysis Details


Analysis Details screen displays all information about the analysis performed with a certain sample.

There are four types of analyses with common parameters as well as specific ones.

 Narrative analyses have a specific field "**Narrative text**" that is considered as a result of the analysis.

 Bottle analyses have an additional field called "**Location**" that does not correspond to the location of the sample, and considered as the result of the analysis as well.

 Single and  multi-component analyses contain components with their own parameters.



\$A_5252 - Spk Amt for Semi-Volatile Organics by SP

Validate Reorder Scratch

Analysis pending
Status

Not valid
Validation state

No value
Validation user

No value
Validation date

APHAM
Analyst

AE46117
Sample ID

B01
VQ ID

051220034
Group ID

No value
Container ID

GENERAL ATTACHMENTS

Analysis name Spk Amt for Semi-Volatile Organics by SP

Assigned instrument IC_ANION

Start date 04/14/2021 05:36 PM

Method reference EPA 525.2 Mod

Assigned analyst APHAM

End date 04/28/2021 05:36 PM

Price 0

Instrument IC_ANION

Due date 05/27/2020 12:00 AM

Comments No value

Preparation date 04/06/2021 05:36 PM

Rlt source No value

Special information form fields

Components

Benzy butyl phthalate From definition Result history Replicates

Bis(2-ethylhexyl)phthalate From definition Result history Replicates


Di-n-octyl phthalate From definition Result history Replicates

Edit

Close

There is a number of buttons in the right-hand side of the screen:

- "**Edit**" button that enables the **Analysis Editing** screen.
- "**Audit Trail**" menu item that allows you to browse all changes made to the analysis on the **Audit Trail** screen.



The detailed information about the results of other samples underwent this analysis is available on the **Result History** screen, activated by the "**Result History**" button next to the corresponding component.

Single and multi-component analyses may have replicates. They are displayed on the **Analysis Replicates** screen accessed by clicking the "**Replicates**" button next to the corresponding component.

The click on the sample **Group ID** name opens the **Log Batch Details** page. If the log batch details do not exist, the **Log Batch Details** page will be open with default values. The new details will be added to the system after saving.

You can perform following actions with the analysis using the button on the left-hand side of the screen:

- Validate (mark as valid)
- Invalidate (mark as invalid)
- Reorder (request reanalysis)
- Scratch (cancel the analysis)

The exact set of available actions depends on the current status of the analysis.

The Attachments related to the Analysis are displayed under the “**Attachments**” tab, with sorting and filtering options. You can manage Attachments on the **Analysis Editing** mode.


Analysis Editing


Analysis Editing mode is enabled when you press the **"Edit"** button on the **Analysis Details** screen or by double-clicking any of the editable fields.

Only users trained to perform an analysis can edit this analysis. This rule applies to the entire site.

In the edit mode, certain fields of the analysis become editable and the **"Edit"** button is replaced with **"Save"** and **"Cancel"** buttons that apply or discard the changes respectively.

Note that every change to the analysis while it is in the **"Validated"** state, changes it to **"Not Valid"**.

Every change can be discarded separately by the **"Revert"** button  next to the updated field.




The screenshot shows the 'Analysis Editing' interface for the analysis titled 'Spk Amt for Drinking Water VOCs - THMs'. The interface is divided into several sections:

- Left Sidebar:** Contains status indicators and validation information. It shows 'Manually entered' status, 'Not valid' validation state, and 'No value' for validation user and date. It also shows 'No value' for validation date and 'No value' for analyst.
- General Section:** Contains fields for 'Analysis name' (Spk Amt for Drinking Water VO), 'Method reference' (EPA 524.2), 'Price' (12), and 'Comments'. It also has fields for 'Assigned instrument', 'Assigned analyst', and 'Instrument'.
- Attachments Section:** Contains fields for 'Start date', 'End date', 'Due date' (05/26/2020 10:40 AM), 'Preparation date', and 'Rit source' (No value).
- Special information form fields:** A section for additional information.
- Components:** A list of components with their definitions and result history/replicates links. The components listed are: Bromodichloromethane, Bromoform, Chloroform, and Dibromochloromethane.

The Analyst can be selected in the user selection window opened by the button next to the field.

The potential analysts are displayed in a paginated list with sorting and filtering options. Only users trained to perform the analysis can be displayed in the list.

To set a user as the analyst, select him in the list and click the **"Apply"** button.



The screenshot shows the 'Analyst' selection window. It features a search bar with a filter icon and a dropdown menu. Below the search bar, there is a table of potential analysts:

Initials	First and last name
AKALINSKY	Arseny Kalinsky
APHAM	Labworks Fuser
LONGNAME1234	FirstNameFirstNameFI LastNameLastNa...

At the bottom of the window, there is a 'Items per page' dropdown set to 10, showing '1-3 of 3 users'. The window also includes 'Cancel' and 'Apply' buttons.

It is possible to add custom components to multi-component analysis, beside the default ones that have "From definition" label on them, by clicking the "Add Component" button and submitting the form with component details in the activated window. The custom components can be renamed or deleted later.

Add Component

Name

Start typing

Result

Secondary result

Qualifier

MDL

Dilution factor

1

PQL

|

Result unit

ug/L

Raw result

Result method

EPA 200.8/6020

Result start date

04/14/2021 06:52 PM

Result end date

04/14/2021 06:52 PM

Cancel

Add

In-line Result Entry

It is possible to enter the analysis results across the application using the in-line result entry function.

Only users trained to perform an analysis can edit this analysis. This rule applies to the entire site.


Note that every change to the analysis while it is in the "**Validated**" state, changes it to "**Not Valid**", although it is not visible in the in-line result entry.

The results for single-component and bottle analyses can be entered directly in the list with the analysis by clicking the result value. To discard or apply the updated value, you have to click the corresponding button appeared to the right.

Result
12 ✖ ✓

The narrative analysis, on the other hand, shows a dialog window when you click on the result value.


You can update the result by changing the narrative text and clicking the "**Save**" button.




 !MICROPS ✕

Narrative text

Cancel Save

The multi-component analyses open a window with a list of components and their results. You can specify the value of each component in the window and apply them by clicking the "**Save**" button.

 \$A_5252 ✕

Component name	Result	
Benzy butyl phthalate	<input type="text"/>	
Bis(2-ethylhexyl)phthalate	<input type="text"/>	
Di-n-octyl phthalate	<input type="text"/>	

Cancel Save

Note that changing analysis result may affect the result of another analysis if it is a part of calculation.

Analysis Result History


You can access the **Result History** screen from analysis details or the enter results table.

The screen contains details about the results of a particular analysis with a list of results of related samples taken in the specified location. Note that only those samples the user has access to are taken into account.

The first part of the screen contains general information about the analysis with the most recent result values displayed in a chart. The data in the chart may be ordered either by **Sample ID** or the **Collected date** field for single- and multi-component analyses. Bottle analyses display up to 20 most recent locations with their corresponding numbers of occurrences in form of a chart. Narrative analyses display a list of results with the number of occurrences of each particular result.

The second part of the screen is a list of affected samples with an option to sort and filter by corresponding fields.


The results are color-coded according to the **Analysis Legend**.



Analysis Replicates

Analysis Replicates screen displays the details about the replicates - a set of repetitions of the analysis intended to test for the deviation of the results.

The **"Edit"** button on the right-hand side of the screen enables the **Analysis Replicates Editing** mode and the **"Audit Trail"** button allows you to browse all changes made to the replicates of a certain analysis on the **Audit Trail** screen.



The replicates can be added in the window that can be accessed by clicking the **"Add replicate"** button in the right-hand side of the **"Replicates"** table.

To be able to add multiple replicates without closing the window, please select the **"Add multiple"** option under the drop-down menu next to the **"Add"** button.

The 'Add Replicate' dialog box has fields for 'Result', 'MDL', 'Dil. factor', and 'Comments'. The 'Dil. factor' field is set to '1'. At the bottom, there are 'Cancel' and 'Add multiple' buttons, with a dropdown arrow next to 'Add multiple'.

Upon addition, the replicate becomes available in the **"Replicates"** table changing the screen to the editing mode automatically.

When the analysis has at least one replicate, the information is also displayed in the chart.

Analysis Replicates Editing

Analysis Replicates Editing mode is enabled when you press the **"Edit"** button on the **Analysis Replicates** screen or by double-clicking any of the editable fields.

In the edit mode, certain fields of the analysis become editable and the **"Edit"** button is replaced with **"Save"** and **"Cancel"** buttons that apply or discard the changes respectively.

The screenshot shows the 'AP_A1 Analysis' window. On the left, there's a sidebar with 'Sample ID: AE45621' and 'Analysis: AP_A1 - AP_A1 Analysis'. Below this, a dropdown menu is set to 'Average replicate captured'. Three radio buttons are present: 'From all available replicates' (selected), 'From selected replicates', and 'From replicates in MDL range'. A 'Default MDL' field is empty. The main area displays 'Analyst: APHAM' and a 'Date' of '04/15/2021 02:45 PM'. A 'Result' field shows '2'. Below this, 'MDL' and 'Dil. factor' are set to '1'. A bar chart on the right shows six replicates with a linear trendline. Below the chart, 'Average' and 'Standard deviations' are calculated. At the bottom, a 'Replicates' table lists six replicates with their respective results, MDLs, dilution factors, and target MDLs. An 'Add replicate' button is in the top right corner.

Replicate	Result	MDL	Dil. factor	Target MDL
1	1	5	1	5
2	2	5	1	5
3	3	5	1	5
4	1	5	1	5
5	2	5	1	5
6	3	5	1	5

The replicates can be added in the window that can be accessed by clicking the **"Add replicate"** button in the right-hand side of the **"Replicates"** table.


To be able to add multiple replicates without closing the window, please select the **"Add multiple"** option under the drop-down menu next to the **"Add"** button.

The **"Default MDL"** field value will be used as the default one for new replicates.

The 'Add Replicate' dialog box is shown. It has fields for 'Result' (set to 1), 'MDL' (empty), 'Dil. factor' (set to 1), and 'Comments' (empty). At the bottom, there are 'Cancel' and 'Add multiple' buttons, with a dropdown arrow next to 'Add multiple'.

Upon addition, the replicate becomes available in the "**Replicates**" table, where you can change the data of each replicate directly anytime. Every change can be discarded separately by the "**Revert**" button next to the updated field in the table.

The drop-down menu to the right of each replicate can be used to set the replicate data as the output data of the analysis.




Use the cross buttons on the left to delete replicates.


When there are multiple replicates, you can also set the output result using the automatic options to the left from the chart:

- Minimum replicate captured - the replicate with minimum result value will be set as the output data.
- Maximum replicate captured - the replicate with maximum result value will be set as the output data.
- Average replicate captured - the average result value will be set as the output data.
- Custom result - the result can be specified manually by typing into the output data fields.

You can also limit the set of involved replicates by taking into account only the replicates in MDL range, or by selecting multiple replicates in the table and checking the "**From selected replicates**" option below the "**Set analysis to**" drop-down list.



The way the replicate data is displayed in the chart can also be configured. To access the configuration window, click the "**Configure**" button under the drop-down menu in the top-right corner of the chart.



You can choose to display up to two trendlines with a number of preset functions, some of which are parametric.

There is an option to display boundary lines corresponding to the minimum, maximum and the average values.

The chart also allows you to configure the error bar display with multiple options.

Configuration

First trendline

Linear

Second trendline

Polynomial

Order

2

Boundary lines

Minimum, average, maximum

Error bar

Fixed

Value

1

Cancel
Apply

When you want to set the intersection of the trendlines as the output result value, you can click a certain point on the chart. The chart can be zoomed in and out using the mouse wheel, and panned by holding the left mouse button.

Enter Results

The **Enter Results** screen provides an option to enter results for multiple analyses at once.

Only users trained to perform an analysis can edit this analysis. This rule applies to the entire site.

The available information on selected samples with corresponding analyses and their results is displayed in the form of a spreadsheet.

Enter Results

Orientation

Highlight

Nothing

Show comments

Template

MTEMP342019

Revert


Save

	WQ ID	Sample location	Collected date	Submitted	Test date	QA/QC batch me...	Properties	<div><div></div>1623C1</div>	<div><div></div>1623C2</div>	<div><div></div>1623C3</div>	<div><div></div>1623CA</div>	<div><div></div>1623CB</div>	<div><div></div>1623CC</div>	<div><div></div></div>
Name								Total Crypto Counted (...)	Total Volume Analyzed...	Cryptosporidium Oocy...	Crypto Total Number D...	Crypto Total Number D...	Crypto Total Number D...	C...
Department								MICROBIOLOGY	MICROBIOLOGY	MICROBIOLOGY	MICROBIOLOGY	MICROBIOLOGY	MICROBIOLOGY	M...
AE46118	B01	Northwood Zon e 2 Booster PS	05/13/2020	05/12/2020 06:51 PM			Result	0.2	1	23				
							Secondary result							
							Qualifier							
							Due date	05/13/2020	05/13/2020	05/13/2020	05/13/2020	05/13/2020	05/13/2020	05/13/2020
							MDL			0.1				
							Start date	04/15/2021	04/15/2021	04/15/2021				
							Start time	09:18 AM	09:18 AM	09:19 AM				
							End date	04/15/2021	04/15/2021	04/15/2021				
							End time	09:18 AM	09:18 AM	09:19 AM				
							Analyst assigned to							
							Analyst	APHAM	APHAM	APHAM				
							CAS number							
							Dil. factor	1	1	1				
							Instrument assigned to							
							PQL							
PAM dil. factor														
PAM reporting limit														
Raw result														
Result method	EPA 1623	EPA 1623												
Result unit	Oocysts	Liter	Oocysts/L											
Validated indicator	No	No	No	No	No	No	No							
AE46134	B01	Northwood Zon e 2 Booster PS	05/13/2020	05/12/2020 06:51 PM			Result							
							Secondary result							
							Qualifier							

You can view the details about the calculation using the same context menu if the analysis participates in the calculation.

The result history can be viewed by clicking on the corresponding icon in the top-right corner of the table. It will be displayed in a sidebar and update accordingly when you select different analysis result cells in the table. It shows information about up to 10 latest samples underwent the analysis.


The "**Statistics**" tab contains only the samples displayed in the table at the moment.



To add one of the displayed analyses to a sample, select the "**Add analysis to sample**" option under the context menu opened by right-clicking at the intersection of the sample and corresponding analysis.

To update analysis comment, select the "**Comment**" option under the context menu opened by right-clicking any of the analysis field values.

You can update any of the displayed field values by typing directly in the spreadsheet. Using the **"Fill single comp results"** dialog window, you can update single-component analyses results in bulk. It can be opened by the corresponding option under the context menu opened by right-clicking on a particular row, column, or cell range. The value entered in this window is applied to all analyses in the selected range.



×

Fill Analysis Results

☐ Fill pending analyses

☒ Fill all analyses

Fill with a value

Start typing

Cancel

Fill

The click on multi-component analysis result value opens a dialog window dedicated to this particular analysis. You can enter the data for any of the components either manually or through **"Fill results for analysis <Analysis code>"** dialog and click the **"Apply"** button to close the window.

Components Show replicates												
	Result	MDL	Secondary result	Qualifier	Raw result	Dil. factor	PQL	Result unit	Analyst	Start date	Start time	End date
1,1,2-Trichloroethane	0.11	0.11				1	0.50		APHAM	04/15/2021	10:00 AM	04/15/2021
1,1-Dichloroethane	0.094	0.094				1	0.50		APHAM	04/15/2021	10:00 AM	04/15/2021
1,1-Dichloroethane	0.10	0.10				1	0.50		APHAM	04/15/2021	10:00 AM	04/15/2021
1,1-Dichloropropene	0.088	0.088				1	0.50		APHAM	04/15/2021	10:00 AM	04/15/2021
1,2,3-Trichlorobenzene	0.11	0.11				1	0.50		APHAM	04/15/2021	10:00 AM	04/15/2021
1,2,3-Trichloropropene	0.087	0.087				1	0.50		APHAM	04/15/2021	10:00 AM	04/15/2021
1,2,4-Trichlorobenzene	0.13	0.13				1	0.50		APHAM	04/15/2021	10:00 AM	04/15/2021
1,2,4-Trimethylbenzene	0.10	0.10				1	0.50		APHAM	04/15/2021	10:00 AM	04/15/2021
1,2-Dibromo-3-chlorobenzene	0.25	0.25				1	0.50		APHAM	04/15/2021	10:00 AM	04/15/2021
1,2-Dibromoethane	0.14	0.14				1	0.50		APHAM	04/15/2021	10:00 AM	04/15/2021
1,2-Dichlorobenzene	0.082	0.082				1	0.50		APHAM	04/15/2021	10:00 AM	04/15/2021
1,2-Dichloroethane	0.035	0.035				1	0.50		APHAM	04/15/2021	10:00 AM	04/15/2021
1,2-Dichloropropane	0.064	0.064				1	0.50		APHAM	04/15/2021	10:00 AM	04/15/2021

You can validate, invalidate, reorder analysis either for a particular analysis or in bulk for an entire sample or analysis code. It can be done by selecting the corresponding option under the context menu opened by right-clicking on a particular row, column, or a cell.

AE12324	AT3_GRAB	MWRP Aeration Tank #3	05/01/2016
---------	----------	--------------------------	------------

Copy
Paste
Open sample details
Fill single comp results for sample AE12324
Validate all completed analyses for sample AE12324
Invalidate all analyses for sample AE12324
Reorder all analyses for sample AE12324
Sheet operations

You can apply the operations to the entire sheet by right-clicking any cell.

AE12326	AT5_GRAB	MWRP Aeration Tank #5	05/01/2016	07/10/2017 06: 13 PM		TS
---------	----------	--------------------------	------------	-------------------------	--	----



Copy
Paste
Open sample details
Fill single comp results for sample AE12326
Validate all completed analyses for sample AE12326
Invalidate all analyses for sample AE12326
Reorder all analyses for sample AE12326
Sheet operations

Fill all single comp results in sheet
Validate all completed analyses in sheet
Invalidate all analyses in sheet
Reorder all analyses in sheet

You can scratch a particular analysis by right-clicking a cell and selecting the “**Scratch**” operation from the context menu.

03:59 PM	
05/02/2016	Copy
03:59 PM	Paste
	Open sample details
	Open analysis details
JSANTY	Open replicates
	Comments
1	Assign
	Validate
	Invalidate
1.0	Reorder
	Scratch
	Sheet operations

You can assign an analyst or an instrument to a particular analysis or to all unassigned samples for a particular analysis code.




All changes made on the Enter Results screen can be saved by clicking the **"Save"** button in the top-right corner of the screen.

To discard the changes, click the **"Revert"** button.


Note that, like with inline result entry, every change made to the analysis while it is in the **"Validated"** state, sets it to **"Not Valid"**, although this is not indicated on **Enter Results** screen.

Edit / Review Samples

Edit / Review Samples screen provides a way to update samples in form of a spreadsheet that may contain one or many samples. The spreadsheet supports MS Office Excel style operations.





Sample ID	Location code	A_TEST	AGENCY	Analyses	Attachment	b_TEST	Billing rate	Check specifications	Collection date	Collection time	Comments	Coordinate
AE17475	AT3_GRAB			Analyses	Attachments		1	<input type="checkbox"/>	07/17/2016	03:17 PM	Comments	Not se
AE17477	AT5_GRAB			Analyses	Attachments		1	<input type="checkbox"/>	07/17/2016	03:19 PM	Comments	Not se
AE17478	AT6_GRAB			Analyses	Attachments		1	<input type="checkbox"/>	07/17/2016	03:20 PM	Comments	Not se

You can navigate to this screen from the **Explorer** or **Work Queue Manager** screens by selecting one or more samples and clicking the  button or by using the context menu.

You can edit / review the samples by selecting the “**Edit / Review**” option under the sample group action menu one of the following screens:

- Log Batch Explorer
- QA/QC Batch Explorer
- Collection Group Explorer

Additionally, you will find the  button on most screens of the Sample Management. You can either click on this icon to navigate to the **Edit / Review Samples** with a template that you used last, or select a particular template by selecting it in the drop-down list next to this button.




The set of fields available for editing depends on the selected template and the field data input works the same way as it does on the **Create Samples** screen.

In case if you want to revert your changes, click the “**Revert changes**” button next to the “**Save**” button.


When you are done editing samples, click the “**Save**” button to apply the changes.

Receive Samples

Receive Samples screen is similar to the **Edit / Review Samples** screen, except that after the sample data is saved, the samples automatically become received.




Sample ID	Location code	A_TEST	AGENCY	Analyses	Attachment	b_TEST	Billing rate	Check specifications	Collection date	Collection time	Comments	Coordi
AE48503	WHATAU		IRWD	Analyses	Attachments		0		04/07/2021	12:00 AM	Comments	Not se
AE48505	WHATAU		IRWD	Analyses	Attachments		0		04/07/2021	12:00 AM	Comments	Not se

You can navigate to this screen from the **Explorer** or **Work Queue Manager** screens by selecting one or more samples and clicking the  button.

You can receive the samples by selecting the “**Receive**” option under the sample group action menu one of the following screens:

- Log Batch Explorer
- QA/QC Batch Explorer
- Collection Group Explorer

Additionally, you will find the  button on most screens of the Sample Management. You can either click on this icon to navigate to the **Receive Samples** with a template that you used last, or select a particular template by selecting it in the drop-down list next to this button.



The “**Receive**” option is available based on the workflow architect rules defined in the State Management desktop application.

When clicking the “**Save and mark received**” button, you will be requested for a reason for the **Audit Trail** entry or **Electronic Signature**. It depends on the workflow architect settings and audit trail configurations.

You can point out to print a custody tracking report or labels after confirming the **Receive** operation.

Electronic Signature

Electronically sign to verify that the data has been changed and verified by you.

Login

APHAM

Password

....

Comments

Start typing


☐ Print custody tracking report

☐ Print labels

Cancel

Sign

If you want to skip the editing process and just mark the sample group as received, you can select the **“Mark received”** option under the sample group action menu on the **Log Batch Explorer**, **QA/QC Batch Explorer**, **Collection Group Explorer** pages.




Validate samples

Validate Samples screen provides an option to validate samples.

Sample ID	Validation availability	Analysis valida...	Violation	Audit trail	QA/QC data	Comments	Special information
AE47202	Ready to validate	Yes	No	Yes	No	Comments	Special info
AE47201	Ready to validate	Yes	No	Yes	No	Comments	Special info
AE47200	Ready to validate	Yes	No	Yes	No	Comments	Special info
AE47199	Ready to validate	Yes	No	Yes	No	Comments	Special info
AE47197	Ready to validate	Yes	No	Yes	No	Comments	Special info

You can navigate to this screen by selecting the “**Validate**” option under the “dot-dot-dot” menu of the **Explorer**, **Work Queue Manager**, **Sample Details** screens. The screen is also available from the **Explorer** and **Work Queue Manager** screens via the context menu for selected samples.



You can navigate to this screen by selecting the “**Validate**” option under the sample group action menu one of the following screens:

- Log Batch Explorer
- QA/QC Batch Explorer
- Collection Group Explorer

The “**Validate**” option is available based on the workflow architect rules defined in the State Management desktop application.

The **Validate Samples** screen consists of two tabs: Samples and Analyses.

The **Samples** tab provides users with the ability to verify sample fields. The set of fields available for viewing and the number of columns to be locked can be configured via the **Settings** dialog available from the gear icon.

If you do not want to validate a sample or it cannot be validated, you can exclude it from the list via the cross icon.

If there are any violations in sample analyses, the **Violation** column will display “Yes”. The click on this value opens the sample details on the Analysis tab so that you can check violation details.

The **Analyses** tab displays the same table as the **Enter Results** screen and provides you with the ability to check and edit analyses details.

LABWORKS

ENTERPRISE LIMS

DASHBOARD

SAMPLE MANAGEMENT

CONFIGURATION

Labworks Fuser

Demo database

Validate Samples

SAMPLES

ANALYSES

Save and validate

Enter Results

Highlight

Nothing

Show comments

Template

Sample_Down

	Properties	AP_A2	B_1623GB	ODOR	PH
AE47183	Result		Absent		
	Status		Manually entered		
	Price		0		
AE47870	Result			Absent	
	Status			Manually entered	
	Price			0	
AE48422	Result	2			7
	Status	Manually entered			Manually entered
	Price	0			0

Unlike the Enter Results screen, you can specify the **Price** value for each analysis.

You can navigate to the details of a particular sample by clicking on corresponding fields within the sample list and the analysis table.

When you are done checking samples and analyses, click the **“Save and Validate”** button to validate the samples. You may be requested for the **Electronic Signature**. It depends on the workflow architect settings.

Scan Samples

The Scan Samples window provides a quick way to find samples using bar code scanner.

Open the Scan Samples window and enter (with or without the bar code scanner) as many Container IDs or Sample IDs as you want to the field in the top-right to add them to the selection list.

Scan Samples						<input type="text" value="Enter sample or container ID"/>	
	Sample ID	Group ID	WQ ID	Collected date	Submitted date	Status	
✖	AE48496	210215181	AQUASPL	02/22/2021 09:59 PM	02/22/2021 09:59 PM	Collected	
✖	AE48493	012221175	WGWTRCA	01/22/2021 12:00 AM	01/22/2021 08:09 PM	Waiting for Valid...	
✖	AE48494	012221175	WGWTRCB	01/22/2021 12:00 AM	01/22/2021 08:09 PM	Waiting for Report	
✖	AE48495	210215181	@0	02/15/2021 04:25 PM	02/15/2021 04:25 PM	Ready to Collect	

When you are done scanning the samples, you can either clear the Sample ID filter on the **Sample Explorer** and apply the selection as is, or add the scanned sample IDs to the current selection using the drop-down menu next to the “Apply selection” button.

Sample ID (selection)	AE48496 +3 x
Sample location	AE48496 AE48493 AE48494 AE48495
Aqua Splash Product	
No value	

The **Sample Explorer** should have these sample IDs displayed in the applied filters' area.

The **Work Queue Manager** simply displays the scanned samples.

Show selected (4)

Apply the **Remove all** button to clear the list of scanned samples in the window.

Audit Trail

The audit trail is a list of changes made to a certain sample, analysis or a replicate, that is displayed in the top right corner of the window.

The elements are displayed can be sorted by clicking the labels in the header of the table. Subsequent clicks reverse the sorting order.

To filter the elements, expand the quick filter bar by clicking the **funnel icon** on the left-hand side of the window, specify the filtering options, and apply them by clicking **check mark icon** in the quick filter bar.

The **cross icon** in the quick filter bar resets the filter.

In order to see the additional events, such as the automatic ones, you can uncheck the "**Display manually entered events only**" checkbox in the top of the window.

Audit Trail Display manually entered events only AE48496						
Changed item	User-friendly changed item	Previous value	New value	Modified date	Changed by	Reason
RU_C24_ECOLVE.coli (Colliert)	C24_ECOLVE.coli (Colliert): Unit	No value	P-A/100mL	03/11/2021 05:51 PM	Labworks Fuser	QC Review
RF_C24_ECOLVE.coli (Colliert)	C24_ECOLVE.coli (Colliert): Result method	No value	SM9223 B	03/11/2021 05:51 PM	Labworks Fuser	QC Review
RD_C24_ECOLVE.coli (Colliert)	C24_ECOLVE.coli (Colliert): Dilution factor	No value	1	03/11/2021 05:51 PM	Labworks Fuser	QC Review
Items per page: 10 1-3 of 3 events						
Close						

Certain actions, such as sample editing or analysis editing may request a reason to be specified.

You can select one of the common reasons by selecting one from the proposed list, or type a custom comment in the window opened after you apply the changes.

Audit Trail

Please specify reason for data change.


Reason

Start typing

CancelSave

Electronic signature

Some operations with objects require confirmation as an electronic signature. If the system is configured to request an electronic signature, the user must enter the correct password for their profile to confirm their action.



The image shows a dialog box titled "Electronic Signature" with a close button (X) in the top right corner. Inside the dialog, there is a message: "Electronically sign to verify that the data has been changed and verified by you." Below this message, there are three input fields. The first is labeled "Login" and has a red border with the text "Please enter your login." below it. The second is labeled "Password" and has a red border with the text "Please enter a password." below it. The third is labeled "Comments" and has a placeholder text "Start typing". At the bottom right of the dialog, there are two buttons: "Cancel" (grey) and "Sign" (orange).

The **Electronic signature** window also requested you a reason for the **Audit Trail** entry.

Electronic signature can be configured by the workflow architect via the desktop State Management application. The default workflow architect template is configured to request the e-signature for the following sample operations:

- Receive
- Validate
- Set sample status to Report completed
- Set sample status to Report rejected
- Set sample status to Invoice completed
- Set sample status to Invoice rejected
- Set sample status to Completed

Analysis Legend







Analysis Types

-  Narrative analysis
-  Single-component analysis
-  Multi-component analysis
-  Bottle analysis
-  Group analysis

Analysis States




-  Validated
-  Not valid

Analysis Statuses


-  Analysis pending
-  Manually entered
-  Uploaded from file
-  Calculated
-  Pre-defined
-  Scratched


Analysis Result Colors

The color of the bar next to the analysis result value depends on the specification violation:


-  The result exceeds specification range
-  The result is exceeding warning range but within the specification range
-  The result is within the warning range

Collection Group Statuses

 Pending (all samples are in Not Ready or Ready to Collect statuses)

 Completed (no samples in Not Ready or Ready to Collect statuses)

 Received (all samples are reviewed)

 In progress (in all other cases)

Checksum Violations

Checksum violation detection module is one of the Labworks Enterprise LIMS security features.

Only the data collected by the Labworks Desktop Application and Labworks Enterprise LIMS is considered as valid. So, when the database is changed directly, the checksum violation detection module notifies the user about the infiltration.

The module covers changes to the most sensitive data:

- Results of analyses
- Audit trail entries

The notification messages are displayed on the **Analysis Details** or **Audit Trail** screens when any modifications to corresponding data have occurred from outside the system.

Loaded data contains one or more check sum violations! The following analysis results have violations:
'Temperature (Field)'



Loaded data contains one or more check sum violations! The following audit trail records have violations:
'MF_TCOL_Total Coliform, Membrane-Filter: Result' modified at '11/17/1997 11:16'

