

---

## **LABWORKS LIMS v6.8**

---

### **Setting up ANALYST in LABWORKS Database**

---

---

**Last Updated:** March 30, 2018


---




## Adding a new ANALYST in LABWORKS Database


### Follow below steps to add analyst in LABWORKS

1. Create new LABWORKS User using System Manager Utility.
2. Launch “Training6.exe” application from “\\LABWORKS\Client” folder.





3. Click on “Load” button. Below “Analyst Selection” window will appear.



Basically, it displays all the Users who are Active or Suspended from USERHEAD table.

4. Select a User from above window to create Analyst.

Personnel Training

Reports Schools Help

AMOL - Amol T

Load Save Delete Exit

Employed On: 03/30/2018  Part Time  Full Time  Lab Analyst?

Load Analysis


Clear All

Anl Code	Analysis name	Expiration	New SOP

Analyses Comments Training Experience License New location list maintenace - Q2 2004

Selected User will be displays on the form.

5. Check “Lab Analyst?” checkbox. To make this user as an Analyst.
6. Click on “Save” button to generate Analyst record for the selected User in “ANALYST” table.



7. Verify "ANALYST" table in database whether analyst is created or not.

SQLQuery1.sql - DTP302.LABWORKS\_6.8\_TEST (sa (73)) - Microsoft SQL Server Management Studio

File Edit View Query Project Debug Tools Window Help

SQLQuery1.sql - DTP...\_6.8\_TEST (sa (73))

```
SELECT * FROM ANALYSTS
```

100 %

	INITIALS	FULLNAME	STATUS	EmpDate	EmpType	Other	hSec	TEAMNAME
1	AMOL	Amol T	1	2018-03-30 00:00:00.000	0		137930	NULL
2	DM	Desk Menu1	1	2018-03-30 00:00:00.000	1		137821	NULL
3	FU		1	2017-07-10 16:30:31.753	0		0	BLUE
4	NU		1	2017-07-10 16:30:06.183	0		0	BLUE
5	RAM	Ram Nikam	1	2018-03-30 00:00:00.000	0		0	NULL
6	SU	Second User	0	2008-05-19 00:00:00.000	0		0	
7	TU		1	2017-07-10 16:29:58.007	0		0	BLUE
8	USR	First User	1	2008-04-29 00:00:00.000	0		0	BLUE

Query executed successfully. DTP302 (12.0 RTM) sa (73) LABWORKS\_6.8\_TEST 00:00:00 8 rows

Output

Ready Ln 2 Col 1 Ch 1 RIS