

The LABWORKS Privilege Matrix, documenting all individual privileges that can be set for each user in version 7.2 of our product.

URVALUE	GENERAL LEVEL	LIMS Administrator	Lab Analyst	Lab Supervisor	Administrative Assistant	IT Administrator	EXPANDED DESCRIPTION
ACCRPSET_APP	ADMIN / USER	X					Allows the user to run the Access Report program
ADD_REPLICATE	ADMIN / USER	X	X	X	X		Allow the user to add new Replicate
AER_RUN	ADMIN	X					Allows the system to run automated email reports
AMAIN_T_APP	ADMIN / USER	X	X	X		X	Allows the user to run the Analysis List Maintenance Program from the LABWORKS Desktop
ANALVALIDATE	ADMIN / USER	X	X	X			Allows the user to validate analyses
ANALYSIS_STATE	ADMIN	X					Allows the user to modify the analysis state in Analysis List Maintenance
ANOBJECTCHANGE	ADMIN / USER	X	X	X			Allows the user to modify or delete an object stored with an analysis
ANOBJECTLOAD	ADMIN / USER	X	X	X			Allows the user to load an object onto an analysis
ANOBJECTVIEW	ADMIN / USER	X	X	X	X		Allows the user to view an object that has been stored on an analysis
APPROVE_PRIV	ADMIN / USER	X		X		X	Run Object Publisher
ARCHIVEV_APP	ADMIN	X					Allows the user to run the Sample Archive Program from the LABWORKS Desktop
ASSIGN_ANALYST	ADMIN / USER	X	X	X			Allows the user to assign an analyst to an analysis

ASSIGN_DUEDATE	ADMIN / USER	X	X	X				Allows the user to assign a due-date on an analysis
ASSIGN_INST	ADMIN / USER	X	X	X				Allows the user to assign an instrument to an analysis
AUDIT_CHKSUMVIEW	ADMIN	X				X		Allows the user to view the Check-sum field within the Audit-trail viewer
AUDIT_FIXCHKSUM	ADMIN	X				X		Allows the user to update the check-sum field on an audit record
AUDIT_PRIV	ADMIN / USER	X				X		Allows the user to access the Audit-trail Viewer
AUDITRUN	ADMIN / USER	X	X	X		X		Allows the user to run the Audit-trail viewer
BACKLOGO_APP	ADMIN / USER	X	X	X				Allows the user to run the Backlog Program from the LABWORKS Desktop
BACKLOGRUN	ADMIN / USER	X	X	X				Allows the user to run a backlog report
BARCODERUN	ADMIN / USER	X	X	X				Allows the user to generate a barcode (within the Sample/Modify Program)
BATCHDELETE	ADMIN	X		X				Delete QA/QC Batches
BATCHSETUP	ADMIN / USER	X	X	X				Create QA/QC Batches
CALCSETUP	ADMIN	X						Allows the user to setup / create calculations in the lims
CAR_MODIFY	ADMIN / USER	X						Allows the user to modify CAPA Reports
CHANGEPASSWORD	ADMIN / USER	X	X	X	X	X		Allows the user to change their password
CINV	ADMIN / USER	X						Allows the user to run the Chemical Inventory program
CINV_APPR_ITEMS	ADMIN / USER	X						Allows user to Approve Item Stock in Chemical Inventory
CINV_CR_ITEMNAME	ADMIN / USER	X						Allows user to create new Item in Chemical Inventory
CINV_CR_VENDOR	ADMIN / USER	X						Allows user to Create Vendor
CINV_MODIFY_ITEM	ADMIN / USER	X						Allows user to edit Chemical Inventory Item Stock

CINV_OVRDEXPDATE	ADMIN / USER	X					Allows user to override expiry date of item in Chemical Inventory.
CINV_PR_APPR	ADMIN / USER	X					Allows user to Approve Purchase Requisition
CINV_RECONCILE	ADMIN / USER	X					Allows user to Physical Reconciliation
CINV_SPUSER	ADMIN / USER	X					Chemical Inventory Super Owner
COA_BYPASS_VIOLS	ADMIN / USER	X					Allows the user to generate a COA that has specification violations
COADELETE	ADMIN	X					Allows the user to delete a COA Customer Base Record
COAPRNO_APP	ADMIN / USER	X					Allows the user to run the COA Program from the LABWORKS Desktop
COARUN	ADMIN / USER	X					Allows the user to run / generate a COA
COASETUP	ADMIN	X					Allows the user to configure the COA program
CONFIGDESKMENU	ADMIN	X				X	Allows the user to configure the LABWORKS Desktop Menu
CONFIGWEBMENU	ADMIN	X				X	Allows the user to configure the LABWORKS Webtop Menu
CONVERTCALDEFS	ADMIN	X					Allows the user to run the ConvertCalc Program (Update calculations from 6.3 to 6.4)
CONVSPEC60TO61	ADMIN	X					Allows the user to run the ConvertSpecs Program (Update specifications from 6.3 to 6.4)
CRPSET_APP	ADMIN	X					Allows the user to run the Crystal Report Setup program from the LABWORKS Desktop
CTRACK_APP	ADMIN / USER	X	X	X	X		Allows the user to run the Custody Tracking program from the LABWORKS Desktop
CUSTODYRUN	ADMIN / USER	X	X	X	X		Allows the user to enter Custody Tracking information into Custody Tracking
DELETE_REPLICATE	ADMIN / USER	X	X	X	X		Allow the user to delete existing Replicate
DESKTOPSETUP	ADMIN	X					Allows the user to register and modify the registration of programs associated with the Labworks Desktop (pre 6 privilege)
DMR_APP	ADMIN / USER	X					Allows the user to run the DMR application from within the LABWORKS Desktop
DMRREAD	ADMIN / USER	X					Allows the user to view DMR information in the DMR Application
DMRRUN	ADMIN / USER	X					Allows the user to create DMR reports information in the DMR Application

DMRSETUP	ADMIN	X					Allows the user to setup DMR information in the DMR Application
DMRWRITE	ADMIN / USER	X					Allows the user to enter DMR information in the DMR Application
DOCREFCREATE	ADMIN	X					Allows the user to create new documents in the Documentation tracking program
DOCREFDELETE	ADMIN / USER	X					Allows user to delete Document Reference
DOCREFO_APP	ADMIN / USER	X	X	X	X		Allows the user to run the Documentation tracking program from the LABWORKS Desktop
DOCREFREVIEW	ADMIN / USER	X					Allows the user to review Documents
DOCREFVIEW	ADMIN / USER	X	X	X	X		Allows the user to view documents referenced in the Documentation tracking program
DV_SHOW_HISTORY	ADMIN / USER	X	X	X	X		Allows the user to Show History from the Dataviewer
ERRORLOGDELETE	ADMIN	X					Allows the user to delete errorlog entries
ERSETUP_APP	ADMIN / USER	X	X	X	X		Allows the user to run the Exception Report setup program from the LABWORKS Desktop
EWRKSHTCOMPLETE	ADMIN / USER	X					Allows the user to complete eWorksheet
EWRKSHTDELETE	ADMIN / USER	X					Allows user to delete eWorksheet
EWRKSHTEDIT	ADMIN / USER	X					Allows the user to Change eWorksheet to In Progress
EWRKSHTSETUP	ADMIN / USER	X					Allows user to setup eWorksheet
EWRKSHTTEMPEXCEL	ADMIN / USER	X					Allows user to Edit eWorksheet excel template
EWRKSHTTEMPSETUP	ADMIN / USER	X					Allows user to setup eWorksheet template
EWRKSHTVALIDATE	ADMIN / USER	X					Allows the user to Validate eWorksheet
EXCEPTIONRUN	ADMIN / USER	X	X	X	X		Allows the user to run Exception Reports from the Exception Report program
EXCPTREPORTSETUP	ADMIN	X					Allows the user to configure the Automatic Exception Reporting program. This program will detect exceptions and report them with a stored Word document, email or sms the exceptions to a user or a group of users

EXCPTRLO_APP	ADMIN / USER	X					Allows the user to run the Automatic Exception Report program from the LABWORKS Desktop
EXPORT_APP	ADMIN / USER	X	X	X	X		Allows the user to run the Export Data program from the LABWORKS Desktop
EXPORTRUN	ADMIN / USER	X	X	X	X		Allows the user to export data from the LIMS using previously configured exporting formats
EXPORTSETUP	ADMIN	X					Allows the user the configure formats for exporting data from the LIMS into some required format. Several options for exporting are available. See the User Manual.
FDETSET_APP	ADMIN / USER	X					Allows the user to run the Field Entry Data Terminal program from the LABWORKS Desktop
FIELDLABELRUN	ADMIN / USER	X					Allows the user to print labels to be used for field sampling.
FORMSETUP	ADMIN	X					Allows the user to create special lab forms using the Forms Generator Program
GROUPTYPECONFIG	ADMIN	X				X	Allows the user to configure Persisted Group information in LWSampleLogin
IMANDC_APP	ADMIN / USER	X	X	X			Allows the user to run the Instrument Maintenance and Calibration program from the LABWORKS Desktop
IMCMODIFY	ADMIN	X	X	X			Allows the user to modify the data within the Instrument and calibration database
INIEDIT	ADMIN	X				X	Allows the user to change the INI (initialization) file
INSTCONV_APP	ADMIN	X					Allows the user to run the Instrument conversion program from the LABWORKS Desktop
INSTLIB_APP	ADMIN	X					Allows the user to run the Instrument Library program from the LABWORKS Desktop
INVOICERUN	ADMIN / USER	X	X	X	X		Allows the user to run / generate invoices using the standard Labworks Invoicing program.
IPPRUN	ADMIN / USER	X					Allows the user to run the Industrial Pretreatment Program. (IPP) This is a US EPA regulatory program that tracks the permits of industry to discharge into a Publicly Owned Water Treatment facility.
IPPSETUP	ADMIN	X					Allows the user to configure the IPP program
LIPPO_APP							Allows the user to run the Labworks Industrial Pretreatment Option
LMAINT_APP	ADMIN / USER	X	X	X			Allows the user to run the Location List Maintenance program from the LABWORKS Desktop
LOCATIONMODIFY	ADMIN	X	X	X			Allows the user to setup the location codes in the Location List Maintenance program.
LOCATIONREAD	ADMIN / USER	X	X	X			Allows the user to view / print location code information from Location List Maintenance
LOGINWOLOCCODE							Allows the user to login samples that do not have a predefined Location Code
LOGPRINTSETUP	ADMIN	X					Allows the user to configure the information that may be printed a sample login
LOGTEMPLATEDEL							Allows the user to delete a Sample Login Template

LOGTEMPLATESETUP	ADMIN	X						Allows the user to configure the single sample, multi-sample and ProcessScheduler sample login screens or templates
LWEXP_INTERNAL	ADMIN/USER	X	X	X	X	X		Labworks Explorer – Users without this priv don't see samples in LWExplorer until they are Validated.
LWEXPCONFIGURE	ADMIN / USER	X	X	X	X			Allows the user to configure LABWORKS Explorer
LWEXPRUN	ADMIN / USER	X	X	X	X	X		Allows the user to run LABWORKS Explorer
LWEXPSHOWRAW	ADMIN / USER	X	X	X				Allows the user to view un-validated data using the read only Labworks Explorer program
LWSQCO_APP	ADMIN / USER	X	X	X	X			Allows the user to run the LWSQC6 program from the LABWORKS Desktop
LWTRENDO_APP	ADMIN / USER	X	X	X	X			Allows the user to run the LWTREND program from the LABWORKS Desktop
LWUPDATES	ADMIN	X					X	Allows the user to apply database updates in the LWMigrationUtility program
LWWEDCEF_APP	ADMIN / USER	X						Allows the user to run the Edit CEF program from the LABWORKS
LWWFRMED_APP	ADMIN / USER	X						Allows the user to run the Form Editor program from the LABWORKS Desktop
LWWINIED_APP	ADMIN	X					X	Allows the user to run the INI Editor program from the LABWORKS Desktop
LWWINV_APP	ADMIN / USER	X	X	X	X			Allows the user to run the Invoice program from the LABWORKS Desktop
LWWMULTI_APP	ADMIN / USER	X	X	X	X			Allows the user to run the Multi-sample Login program from the LABWORKS Desktop
LWWSINGL_APP	ADMIN / USER	X	X	X	X			Allows the user to run the Single-sample login program from the LABWORKS Desktop
LWWSYSTM_APP	ADMIN	X					X	Allows the user to run the System Manager program from the LABWORKS Desktop
MAILLISTSETUP	ADMIN	X	X	X	X			Allows the user to setup mail-list records in SampMaillist
MAKETEMPLATE	ADMIN	X		X			X	Allows user to create new templates under Demographic Data Entry
MCXFERO_APP	ADMIN / USER	X	X	X				Allows the user to run the multi-component transfer program from the LABWORKS Desktop
METHODREFEDIT	ADMIN / USER	X						Allows user to add custome Method Reference.
MIGRATE_DATA	ADMIN / USER	X						Allows the user to Run migration tool application

MIXSETUP	ADMIN	X					Allows the user to setup special files for calibration/spike mixture files
MODIFY_REPLICATE	ADMIN / USER	X	X	X	X		Allow the user to change Replicate policy
MOR_RPT_APP	ADMIN / USER	X	X	X	X		Allows the user to run the MOR Report program from the LABWORKS Desktop
MORRUN	ADMIN / USER	X	X	X	X		Allows the user to run / generate MOR Reports from the MOR Report program
MORSETUP	ADMIN	X	X	X			Allows the user to configure the MOR (Exel based reporting tool)
MORSTYLE_APP	ADMIN	X	X	X			Allows the user to run the MOR Setup program from the LABWORKS Desktop
MSAMPLES_APP	ADMIN / USER	X	X	X			Allows the user to run the Modify/Delete Samples program from the LABWORKS Desktop
NGIICREATEMODIFY	ADMIN	X					Allows the user to modify/create NGII Instrument Interfaces
NULL_REPL_VALUE	ADMIN / USER	X					Allows user to post Null Replicate result value
PAMRUN	ADMIN	X					Allows the user to use the Project Analyte Management (PAM) program to specify reporting limits on samples
PAMSETUP	ADMIN	X					Allows the user to configure the PAM program
PANLMGTO_APP							Allows the user to run the Project Analyte Management Program
PDETSETUP	ADMIN / USER	X					Allows the user to use the Portable Data Entry program (PDET)
PG_PRIV							Allows the user to access the Persisted Groups in LWSampleLogin
PICKLISTSETUP	ADMIN	X					Allows the user to create picklists for use in certain Labworks programs and forms
PICKLSTO_APP	ADMIN / USER	X	X	X	X		Allows the user to run the Picklist Maintenance program from the LABWORKS Desktop
POST_REPL_VIOL	ADMIN / USER	X					Allows user to post Replicate results with violation
PRELOG_MODIFY							Allows the user to modify pre-login group information
PREPORT_APP	ADMIN / USER	X	X	X	X		Allows the user to run the Progress Report program from the LABWORKS Desktop
PRICELISTDELETE	ADMIN	X					Allows the user to delete a price list in the Location List Maintenance program
PRICELISTSETUP	ADMIN	X					Allows the user to create special prices lists for use in the invoice program. Linked to special clients and customers so that discounts and or special billing can be managed.
PRINT_FROM_MOD	ADMIN / USER	X	X	X			Allows the user to print routing and work sheets from the modify and delete program

PROGRESSRUN	ADMIN / USER	X	X	X	X			Allows the user to run a Progress Report from the Progress Report program
PROGRESSETUP	ADMIN	X	X	X				Allows the user to configure the Progress Report templates for customized viewing of the Labworks data
PROJECT_MODIFY								Allows the user to modify Projects
PROJECT_DELETE								Allows the user to delete Projects
PSEDITSCHEDULE	ADMIN	X		X				Allows the user to configure the Process Scheduler scheduling tool to create an automated sampling plan.
PSLABRECEIPT	ADMIN / USER	X	X	X	X			Allows the user to acknowledge Laboratory Receipt of samples in Process Scheduler
PSPRELABDELETE	ADMIN / USER	X	X	X				Allows the user to delete samples in Process Scheduler that do not have a LABWORKS Sample number assigned
PSPRELABMODIFY	ADMIN / USER	X	X	X	X			Allows the user to modify samples in Process Scheduler before a LABWORKS Sample number is assigned
PSPRODUCTOK	ADMIN / USER	X	X	X	X			Allows the user to acknowledge Production OK of samples in Process Scheduler
PSREORDER	ADMIN / USER	X						Allows the user to re-order samples within Process Scheduler
PSRESET_COLUMN	ADMIN / USER	X						Allows the user to reset the column configuration within Process Scheduler
PSSAMPLELOGIN	ADMIN / USER	X	X	X	X			Allows the user to login samples within Process Scheduler
PSSHELFLIFEDEL	ADMIN	X						Allows the user to delete Shelf-life groups within Process Scheduler
PSSHELFLIFEEDT	ADMIN	X						Allows the user to edit Shelf-life groups within Process Scheduler
PSSHELFLIFEMAK	ADMIN / USER	X						Allows the user to make Shelf-life groups within Process Scheduler
PTMODIFY	ADMIN	X						Allows the user to modify the data within the Personnel Training database
QABATCHE_APP	ADMIN / USER	X	X	X				Allows the user to run the QA Batch Builder within the LABWORKS Desktop
QAQCREPO_APP	ADMIN / USER	X	X	X				Allows the user to run the QAQC Report program from the LABWORKS Desktop
QATEMP_APP	ADMIN	X	X	X				Allows the user to run the QA Template Setup program from the LABWORKS Desktop
QATEMPSETUP	ADMIN	X	X	X				Allows the user to setup QA Templates from the QA Template Setup program
QE_DELETEQUERY	ADMIN	X						Allows the user to delete previously saved Query configured using the Query Engine

QE_EDITDATAMODEL	ADMIN	X					Allows the user to modify an existing Data Model (collection of tables containing the data being searched and the relationships required to allow that data to be retrieved accurately). This requires a person to understand the data structure and to have some knowledge of relational databases.
QE_EDITQUERY	ADMIN	X					Allows the user to modify an existing Query to change the output of the query
QE_NEWDATAMODEL	ADMIN	X					Allows the user to create a new data model by selecting the tables required to populate the queries that will be designed to extract data from the LIMS. This requires knowledge of the data structure and a basic understanding of relational databases.
QE_NEWQUERY	ADMIN	X					Allows the user to create a new query based upon an existing Data Model.
QE_RUNQUERY	ADMIN / USER	X					Allows the user to run an existing query from the Query Engine
QENGINE_APP	ADMIN / USER	X					Allows the user to run the Query Engine program from the LABWORKS Desktop
REENTRY_APP	ADMIN / USER	X	X	X	X		Allows the user to run the Result Entry program from the LABWORKS Desktop
REPORTRUN	ADMIN / USER	X	X	X	X		Allows the user to generate final reports using Labworks standard reporting tools
REPORTSETUP	ADMIN	X					LEGACY PRIVILEGE THAT DOES NOT USE CURRENT REPORTING TOOLS
RESET_CONFIG	ADMIN	X					Allows the user to Reset configuration of all users
RESULTCHANGE	ADMIN / USER	X	X	X			Allows the user to modify a result value that has been saved into the database.
RESULTENTER	ADMIN / USER	X	X	X			Allows the user to load data, import data and calculate results. Required for almost all Labworks users who can add results to a sample
RESULTMERGE	ADMIN / USER	X					Allows user to merge previously posted result
RLTEMPLATESETUP	ADMIN	X					Allows the user to create, modify and configure new or existing result entry templates to meet a users need for a custom data entry screen
RPT_TMPLT_STATE							Allows the user to update the Report Template State
RSETUP_APP	ADMIN / USER	X	X	X			Allows the user to run the Result Entry Template program from the LABWORKS Desktop
SAMPARCHIVE	ADMIN	X					Allows the user to archive data using the Labworks data archival program
SAMPCOLLECTGRP	ADMIN / USER	X					Allows user to create or receive sample collection groups
SAMPDELETE	ADMIN	X					Allows the user to delete samples from the database. Once deleted, there is no trail for that sample. The audit trail ends with the sample being "deleted"

SAMPINACTIVATE	ADMIN	X	X	X	X			Allows the user to set a sample to an inactive status
SAMPLOCKOVERRIDE	ADMIN	X		X				Allows the user to modify stored sample and result information if the LOCKVALIDSAMPS System Manager key is turned on.
SAMPLOGIN	ADMIN / USER	X	X	X	X			Allows the user to create samples in the LIMS. This is not the same as adding a sample at "not ready to collect" or at "ready to collect" status in Process Scheduler
SAMPMODIFY	ADMIN / USER	X	X	X				Allows the user to modify a sample after it has been logged into the system
SAMPRESTORE	ADMIN	X						Allows the user to restore a sample to the active database using the Labworks Restore program
SAMPVALIDATE	ADMIN / USER	X		X				Allows the user to approve the data prior to making it available for reporting or for a Certificate of Analysis
SAVEDESKTOP								LEGACY PRIVILEGE THAT IS NOT CURRENTLY USED
SAWSCREATE								CUSTOM Program for a customer, no general applicability
SAWSSCHVIEW								CUSTOM Program for a customer, no general applicability
SCANDAEMON_APP	ADMIN	X				X		Allows the user to run the ScanDaemon program from the LABWORKS Desktop
SDGROUPSETUP	ADMIN	X				X		Allows the user to create and modify ScanDaemon Groups
SEND_EMAIL	ADMIN	X				X		Allows the user to send an email through the LABWORKS Desktop
SEND_MESSAGE	ADMIN / USER	X	X	X	X	X		Allows the user to send a message through the LABWORKS Desktop
SET_EMAIL_SERVER	ADMIN	X				X		Allows the user to configure email settings in the LABWORKS Desktop
SI_SETUP_APP	ADMIN	X						Allows the user to run the Special Information Form Designer from the LABWORKS Desktop
SMLENTRYMODIFY								Allows the user to modify records in the Sample Mail List Entry
SPECIALINFOSETUP	ADMIN	X						Allows the user to configure a Special Information Form (SIFO) which can be used at the sample level and at the analysis code level
SPECSETUP	ADMIN	X						Allows the user to set the specifications for location code analysis codes. Location code only to be sample specific or at the Primary Analysis code level to be a global specification for all samples using that analysis code
SQ_CONFIGURE	ADMIN / USER	X	X	X				Allows the user to configure Stepped Queries
SQ_RUN	ADMIN / USER	X						Allows user to run stepped queries
SQ_SAVE_PRIVATE	ADMIN / USER	X	X	X				Allows the user to save 'Stepped Queries' to their Private profile

SQ_SAVE_PUBLIC	ADMIN / USER	X	X	X				Allows the user to save 'Stepped Queries' for use by any LABWORKS User
SQCRUN	ADMIN / USER	X	X	X	X			Allows the user to generate SQC Charts or use the SQC program as configured.
SQCSETUP	ADMIN	X		X				Allows the user to configure the SQC Charting program for all normal Labworks SQC configuration attributes.
STATUSRUN	ADMIN	X						Allows the user to generate the system status report. A very general summary of the LIMS sample volume
STDRESULTSETUP	ADMIN	X						Allows the user to modify the system result abbreviaions list
SYS_CONFIGMENU	ADMIN	X					X	Allows the user to modify the System Menu's in the LABWORKS Desktop
TESTMODIFY	ADMIN	X	X	X				Allows the user to modify the definition of analysis codes that are previously configured in the database
TESTREAD	ADMIN / USER	X	X	X	X			Allows the user to view or print the details of an analysis code
TRAINING_APP	ADMIN / USER	X						Allows the user to run the Personell Training program from the LABWORKS Desktop
TRENDRUN	ADMIN / USER	X	X	X	X			Allows the user to activate and view a Trend Chart display of data from the LIMS in real time
TRENDSETUP	ADMIN	X	X	X				Allows the user to setup a real time chart display of a sample point and it's result data to be viewed on a continuous basis.
TTEDIT	ADMIN	X						Allows the user to run the Tool Tip Editor (TTE). This will enable the changing of the help tips that appear as the mouse will cross over an active command of the program.
USE_CONTAINER	ADMIN / USER	X						Allows the user to manage sample container
USR_CONFIGMENU	ADMIN / USER	X						Allows the user to modify their personal menu in the LABWORKS Desktop
VALIDATEDATA								Allows the user to validate sample data
VIEW_SETUP	ADMIN	X		X		X		Allows the user to select the view template configuration in DataViewer and LWExplorer
VSAMPS_APP	ADMIN / USER	X	X	X				Allows the user the ability to run the Sample and Analysis Validation program from the LABWORKS Desktop
WORKFLOW_PRIV	ADMIN	X					X	Allows the user the ability to run the Workflow Architect in the LABWORKS Desktop
WRDSETUP_APP	ADMIN	X						Allows the user to run the Word Report Setup program from the LABWORKS Desktop
WSBATCHSETUP	ADMIN / USER	X						Allows the user to select analysis codes for batched worksheet generation

XLALCO_APP	ADMIN / USER	X	X	X				Allows the user to run the Excel Calculation Setup program from the LABWORKS Desktop
XLALCRUN	ADMIN / USER	X	X	X				Allows the user to enter results that require an Excel calculation. All results entering persons will require this privilege
XLALCSETUP	ADMIN	X		X				Allows the user to configure the Excel calculation engine to perform very special and complex calculations
XREFO_APP	ADMIN / USER	X	X	X	X			Allows the user to run the Cross Reference Search program from the LABWORKS Desktop
XREFRUN	ADMIN / USER	X	X	X	X			Allows the user to execute a previously configured and saved Cross Reference Search.
XREFSETUP	ADMIN	X		X				Allows the user to configure new Cross Reference Searches used to extract data from the LIMS